

# Green 16™

## User Manual

Model : PHT-65LHS  
Version : 1.44

---

• ENGLISH



Full Version

**vatech**



---

## **i**nnovation **i**nside

"i" stands for 'innovation', one of the core values of VATECH, which aims to expand accessibility of medical solutions to more people.



## Notice

Thank you for purchasing the **Green16 (PHT-65LHS)** extra-oral imaging system.

The primary mission of **Green CT** is to aid dental professionals in providing excellent care in a safe environment that promotes healing. However, **Green CT** is not a product name or a brand name.

**Green16 (PHT-65LHS)** is an advanced digital diagnostic system that incorporates PANO, CEPH (Optional), CBCT and 3D MODEL Scan imaging capabilities into a single system.

This manual describes how to operate the **Green16 (PHT-65LHS)** system. It is recommended that you thoroughly familiarize yourself with this manual to make the most effective use of this equipment.

Observe all cautions, safety messages and warnings which appear in this manual.

Due to constant technological improvement, the manual may not contain the most updated information and is subject to change without prior notice to the persons concerned. For further information not covered in this manual, please contact us at:

**VATECH Co., Ltd.**

**Phone: +82-1588-9510**

**E-mail: [gcs@vatech.co.kr](mailto:gcs@vatech.co.kr)**

This document is originally written in English.

**Green16 (PHT-65LHS)** is referred to as “equipment” in this manual.

**Manual Name: Green16 (PHT-65LHS) User Manual**

**Version: 1.44**

**Publication Date: 2020-02**

**Document Number: VDH-UM-068(16)**

# Table of Contents

<b>Notice</b>	<b>v</b>
<b>Table of Contents</b>	<b>vi</b>
<b>1. Introduction</b>	<b>1</b>
1.1 Overview .....	1
1.2 Indications for Use .....	1
1.3 Intended Purposes .....	1
1.4 Intended User Profiles .....	2
<b>2. General Information</b>	<b>3</b>
2.1 Manufacturer's Liability .....	3
2.2 Owner and Operator's Obligations .....	3
2.3 Conventions in this Manual .....	4
2.4 Marks and Symbols .....	5
<b>3. Warnings and Precautions</b>	<b>9</b>
3.1 General Safety Guidelines .....	9
3.2 Electricity-related Safety Precautions .....	13
3.3 Radiation Safety .....	15
3.4 Warnings .....	16
<b>4. Imaging System Overview</b>	<b>19</b>
4.1 System Components .....	19
4.2 Features .....	19
4.3 Imaging System Options .....	19
4.4 Standards and Regulations .....	20
4.5 Operating Principles .....	20
4.6 Imaging System Configuration .....	21
4.7 Equipment Overview .....	22
<b>5. Imaging Software Overview</b>	<b>37</b>
5.1 PC Specifications (Recommended) .....	37
5.2 EzDent-i / EasyDent .....	38

5.3	Console Software .....	39
<b>6.</b>	<b>Getting Started</b>	<b>43</b>
6.1	Turning on the Equipment .....	43
6.2	Running the Image Viewer (EzDent-i / EasyDent).....	44
6.3	Initiating the Console Software .....	49
<b>7.</b>	<b>Acquiring PANO Images</b>	<b>51</b>
7.1	PANO Imaging Program Overview .....	51
7.2	Configuring Exposure Parameters .....	56
7.3	Patient Positioning .....	63
7.4	X-ray Exposure .....	78
7.5	Finishing the Scan .....	79
7.6	Checking the Captured Images .....	79
<b>8.</b>	<b>Acquiring CEPH Images (Optional)</b>	<b>81</b>
8.1	CEPH Imaging Program Overview .....	81
8.2	Configuring Exposure Parameters .....	83
8.3	Patient Positioning .....	87
8.4	X-ray Exposure .....	97
8.5	Finishing the Scan .....	98
8.6	Checking the Captured Images .....	98
<b>9.</b>	<b>Acquiring CBCT Images</b>	<b>99</b>
9.1	CBCT Imaging Program Overview .....	99
9.2	Configuring Exposure Parameters .....	103
9.3	Obtaining Double Scan Image (optional) .....	109
9.4	Patient Positioning .....	115
9.5	X-ray Exposure .....	122
9.6	Finishing the Scan .....	123
9.7	Checking the Captured Images .....	123
<b>10.</b>	<b>Acquiring 3D MODEL Scan Images</b>	<b>125</b>
10.1	3D MODEL Scan Imaging Program Overview .....	125
10.2	Configuring Exposure Parameters .....	126
10.3	MODEL Positioning .....	129

	10.4	X-ray Exposure .....	130
	10.5	Checking the Captured Images .....	131
<b>11.</b>		<b>Troubleshooting</b>	<b>133</b>
	11.1	Troubleshooting .....	133
	11.2	Error Codes .....	134
<b>12.</b>		<b>Cleaning and Maintenance</b>	<b>153</b>
	12.1	Cleaning .....	153
	12.2	Maintenance .....	154
<b>13.</b>		<b>Disposing of the Equipment</b>	<b>157</b>
<b>14.</b>		<b>Technical Specifications</b>	<b>159</b>
	14.1	Mechanical Specifications .....	159
	14.2	Technical Specifications.....	161
	14.3	Electrical Specifications .....	166
	14.4	Environmental Specifications .....	167
<b>15.</b>		<b>Appendices</b>	<b>169</b>
	15.1	Recommended X-ray Exposure Tables .....	169
	15.2	X-ray Dose Data.....	179
	15.3	Electromagnetic Compatibility (EMC) Information .....	188
	15.4	Hand-wrist Image Evaluation References.....	193
	15.5	Acquiring Images for Pediatric Dental Patients .....	197
	15.6	Abbreviations .....	208

# 1. Introduction

## 1.1 Overview

**Green16 (PHT-65LHS)** is an advanced 4-in-1 digital X-ray imaging system that incorporates PANO, CEPH (Optional), CBCT and 3D MODEL Scan imaging capabilities into a single system.

**Green16 (PHT-65LHS)**, a digital radiographic imaging system, acquires and processes multi-FOV diagnostic images for dentists. Designed explicitly for dental radiography, **Green16 (PHT-65LHS)** is a complete digital X-ray system equipped with imaging viewers, an X-ray generator, and a dedicated SSXI detector.

The digital CBCT system is based on a CMOS digital X-ray detector. The CMOS CT detector is used to capture 3D radiographic images of the head, neck, oral surgery, implant, and orthodontic treatment.

**Green16 (PHT-65LHS)** can also acquire 2D diagnostic image data in conventional panoramic and cephalometric modes.

## 1.2 Indications for Use

**Green16 (PHT-65LHS)** is intended to produce panoramic, cephalometric or 3D digital x-ray images. It provides diagnostic details of the dento-maxillofacial, ENT, sinus and TMJ for adult and pediatric patients. The system also utilizes carpal images for orthodontic treatment. The device is to be operated by healthcare professionals.

## 1.3 Intended Purposes

- Determination of the extent of lesions, tumors, cysts, etc., which cannot be adequately visualized on plain films
- Diagnosis of foreign bodies or displaced roots involving the maxillary sinus
- Diagnosis of bone diseases, cysts, etc., affecting the temporomandibular joints
- Identifying the relationship of the inferior dental canal to a tooth/lesion that is to be removed
- Assessment of fractures on maxilla, mandible, condylar neck and fractures of teeth where plain film imaging is equivocal
- Visualization of the 3D anatomy of the alveolar clefts
- Diagnosis of un-erupted teeth impacted teeth and odontomas
- Diagnosis of root resorption of teeth
- Assessment of cleft palate
- Instant diagnosis of CRS (chronic rhinosinusitis)
- Examinations of the airways for measuring the volume and dimensions of air passages

- Reconstruction of position, malformations, and fractures of maxilla & mandible bones, nasal bone and paranasal sinuses as 3D pictures for operational planning and patient education
- Planning any surgery where 3D analysis of the jaw is required
- Storing Plaster Casts in 3D data

## 1.4 Intended User Profiles

Considerations	Requirement Description
Education	<ul style="list-style-type: none"> <li>▪ Licensed dentists or dental hygienists, radiologists and graduates of relevant bachelor's degree (national qualifications)</li> </ul>
Knowledge	<ul style="list-style-type: none"> <li>▪ Understanding the treatment and diagnosis of dental disease</li> <li>▪ Understanding the terms and guidance of hardware and software of a diagnostic medical radiation device and recognizing device connection, installation, operating conditions</li> </ul>
Language understanding	<ul style="list-style-type: none"> <li>▪ Understanding how to use manuals (English / Korean) or</li> <li>▪ Understanding other language provided</li> </ul>
Experience	<ul style="list-style-type: none"> <li>▪ Understanding the objectives and effects of the diagnosis and treatment of dental disease using diagnostic medical radiation devices</li> <li>▪ Understanding of the normal operation of diagnostic medical radiation equipment</li> <li>▪ Understanding the contents of the <b>User Manual</b></li> </ul>

**IMPORTANT**

Qualified personnel should use the dental X-ray CT (dentists, dental hygienists or radiologists only).

## 2. General Information

### 2.1 Manufacturer's Liability

The manufacturers and retailers of this X-ray equipment assume responsibility for the safe and healthy operation of this product only when:

- A **VATECH**-authorized technician has installed the equipment.
- The equipment has been installed by all the cautions and conditions required for installation.
- The good **VATECH** approved equipment and components have always been used.
- A **VATECH**-authorized agent has performed all maintenance and repairs.
- The **User Manual** has typically used the equipment.
- The equipment damage or malfunction is not the result of an error on the part of the owner or the operator.

### 2.2 Owner and Operator's Obligations

- The owner of this equipment shall perform constancy tests at regular intervals to ensure patient and operator safety. Local X-ray safety regulations must perform these tests.
- The owner of this equipment shall perform regular inspection and maintenance of the mechanical and electrical components in this equipment to ensure safe and consistent operation (IEC 60601-1).

The owner of this equipment shall ensure inspection and cleaning work is performed by the maintenance schedule outlined in **Chapter 12: Cleaning and Maintenance**.

## 2.3 Conventions in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions accompanied.

To prevent physical injury and damage to the equipment, please observe all warnings and safety information included in this document.

	<b>WARNING</b>	Indicates information that should be followed with the utmost care. Failure to comply with a warning may result in severe damage to the equipment or physical injury to the operator and/or patient.
	<b>CAUTION</b>	Indicates a situation that demands prompt and careful action, a specific remedy, or emergency attention.
	<b>IMPORTANT</b>	Indicates a situation or action that could potentially cause problems to the equipment and/or its operation.
	<b>NOTE</b>	Emphasizes important information or provides useful tips and hints.
	<b>RADIATION</b>	Indicates a possible danger from exposure to radiation.
	<b>SINGLE USE</b>	Indicates a component that must be replaced for each new patient.
	<b>ESD susceptibility</b>	Indicates that an item is susceptible to damage from electrostatic discharges.

## 2.4 Marks and Symbols

Symbols	Description	Location
	Dangerous voltage	Power board / Inverter board / Monoblock
	Protective earth (Ground)	Column
	Off (power: disconnected to the <b>Main Power Switch</b> )	Main Power Switch
	On (power: connected to the <b>Main Power Switch</b> )	Main Power Switch
	Alternate current	Label
	Type B Applied Equipment (IEC 60601-1: Degree of protection against leakage current and electric shock: Class 1 equipment)	Label
	Radiation hazard	Label
	Indicates the authorized representative in the European Community.	Label
	The CE symbol indicates that this product complies with the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.	Label
	CSA mark No.266436	Label
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.	Label

## 2. General Information

Symbols	Description	Location
	Addresses where the equipment was manufactured.	Label
	Indicates that electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately.	Label
	Warns ESD hazard.	MCU board / Board package
	Indicates that this equipment is classified as a CLASS 1 LASER PRODUCT by IEC 60825-1 ED.2 regulations.	Label
	Indicates that the user needs to refer to the <b>Instruction Manual</b> .	Label
	Indicates the date of manufacture.	Label
	Indicates the manufacturer's serial number so that the specific equipment can be identified.	Label

### 2.4.1 Label Locations

The label is attached on the right side of the equipment, and it consists of 6 parts as below.

**vatech** ①

Product : Computed Tomography X-ray System  
 Model : PHT-65LHS  
 Power Input : 100-240 V~, 50/60 Hz, 2.0 kVA  
 This X-ray equipment complies with 21 CFR Subchapter J

Mode of operation : Continuous operation with intermittent loading – Needs waiting time (at least 60 times the exposure time) before the next exposure begins  
 Mode de fonctionnement : Fonctionnement continu avec chargement intermittent - A besoin de temps d'attente (au moins 60 fois le temps d'exposition) avant le début de l'exposition suivante

WARNING : X-ray unit may be dangerous to PATIENT and OPERATOR unless safe exposure factors, operating instructions and maintenance schedules are observed.  
 AVERTISSEMENT : Cet équipement à rayons X peut être dangereux pour les PATIENTS et les OPERATEURS si les facteurs d'exposition sécuritaires, les instructions de fonctionnement et les programmes de maintenance ne sont pas respectés.

13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, KOREA  
 VATECH Co., Ltd. Website : www.vatech.co.kr  
 Suite 3, Ground Floor, Chancery House, St. Nicholas Way, Sutton, SM1 1JB UK  
 Vatech Dental Manufacturing Ltd.  
 MADE IN KOREA

② **CAUTION (PRUDENCE)**  
 X-RAY / ATTENTION :  
 X-RAY ON WHEN EQUIPMENT IN OPERATION  
 X-RAY / ATTENTION :  
 X-RAY ACTIVE LORSQUE L'EQUIPEMENT EST EN FONCTIONNEMENT.

③ **X-RAY GENERATOR**  
 Model : DG-07E22T2  
 X-ray Tube Type : D-0525B / Canon  
 Focal Spot : 0.5 mm x 0.5 mm (IEC60336)  
 Output : Max. 99 kV, Max. 16 mA  
 Inherent Filtration : 0.8 mm Al / 50 kV  
 Added Filtration : 1.5 mm Al  
 Total Filtration : Min. 2.5 mm Al

④ **CLASS 1 LASER PRODUCT**  
 The laser diode, Class 1 complies with 21 CFR 1040.10 and 1040.11 except for deviations pursuant to laser notice No. 50, dated June 24, 2007, classified to IEC 60825-1 ED.2

⑤ SN  
 Total weight : Max. XXX kg (Without Ceph: Max. XXX kg);

⑥ For Canada only (with CEPH)

⑥ For Canada only (without CEPH)

**CAUTION : X-RAYS**  
**ATTENTION : RAYONS X**

Device Identifier  
Green 16 SC

**CAUTION : X-RAYS**  
**ATTENTION : RAYONS X**

Device Identifier  
Green 16 SP

## 2. General Information

No.	Item
1	<b>Green16 (PHT-65LHS) Main Label</b>
2	<b>CAUTION Label</b> - X-ray / Attention: X-ray on when equipment in operation.
3	<b>X-RAY GENERATOR Label</b> : 1.6 kW Generator
4	<b>CLASS 1 LASER PRODUCT Label</b>
5	<b>Manufacturer Label</b> : The date of manufacture / Serial Number / Weight of the equipment
6	<b>Device Identifier Label</b> (For Canada only)  For Canada, the Model is distinguished by Green16 SC and Green16 SP. <ul style="list-style-type: none"><li>● Green16 SC: CEPH included</li><li>● Green16 SP: CEPH not included</li></ul>

## 3. Warnings and Precautions



Be sure to strictly observe all warnings and safety instructions included in this manual.



This x-ray unit may be dangerous to patients and operators unless safe exposure factors, operating instructions, and maintenance schedules are observed.

### 3.1 General Safety Guidelines

#### Operator qualifications

This equipment may only be operated by personnel adequately trained in its operation.

- To operate this equipment, the operator must:
  - Have read and understood the **User Manual**.
  - Be familiar with the fundamental structure and functions of this equipment.
  - Be able to recognize the intermittent operation of this equipment and implement appropriate measures to remedy such irregularities.

#### General safety precautions

- Follow the instructions specified in this manual to ensure the safety of both the patient and the operator.
- The operator must always maintain vocal/visual contact with the patient during imaging.
- Do not open or remove the cover panels on this equipment. Always have a trained and authorized service technician to carry out inspection and maintenance of this equipment.
- Do not place any heavy objects on this equipment at any time.
- Do not place any objects within this equipment's field of operation. It may cause property damage.
- Do not push or pull the equipment. Overbalances of the equipment may cause the risk of physical injuries or property damage.
- The operator must instruct the patient to remain still until the equipment arm has stopped moving, and the reset motion is completed.
- Observe all local fire regulations. Always keep a fire extinguisher near the equipment.

### 3. Warnings and Precautions

- The operator of this equipment must be familiar with this equipment's emergency protocols.
- Ensure that this equipment is kept away from water, moisture, or foreign substances always.
- If this product is exposed to water, moisture, or a foreign substance, immediately turn off the main power of the equipment and contact your **VATECH** technical support representative.
- If there are signs of oil leakage, immediately cease all operations of this equipment and contact your **VATECH** technical support representative.
- External equipment intended for connection to signal input, signal output or other connectors, shall comply with relevant IEC Standard (e.g., IEC 60950 for IT equipment and IEC 60601-1series for medical electrical equipment).
- Also, all such combination-system-shall comply with the standard IEC 60601-1, and IEC 60601-1-1 harmonized national standard or the combination. If, in doubt, contact a qualified technician or your local representative.
- Any person or organization that installs an external door interlock switch is responsible for ensuring that it has a radiation indicator or equivalent alarm system to show the state of the current.

#### Ventilation

- Do not close the equipment's ventilation slots in any case. The obstruction of ventilation could result in the equipment overheating due to a lack of air circulation.
- Do not spray any liquid or disinfectant on this equipment. The penetration of these substances may damage the electrical and mechanical components inside. Use a soft cloth to disinfect the ventilation slots.
- Always leave enough space around the PC to allow for proper ventilation.

#### Hygiene



Always disconnect the equipment from the power outlet when disinfecting the surfaces of the equipment.

Never expose this equipment to liquids, mists or sprays. Exposing this equipment to liquids may cause an electric shock or otherwise damage the system.

Do not use spray cleaners on the equipment, as this could cause a fire.

- All movable patient support components (the Bite, the Chinrest, the Temple Supports, and the Ear Rods) can be cleaned using alcohol-based cleaning solutions.
- Clean the Support Handles by using alcohol-based cleaning solutions before taking photos of the next patient.
- Other surfaces of the equipment, including the Touch Screen, can be cleaned using a soft cloth dampened with a mild cleaning solution.
- New hygiene cover must be provided for each new patient to prevent the transmission of communicable diseases.



Do not use aerosol or spray cleaning agents directly on the surface of the equipment.

### **Condensation**

- Extreme fluctuation in temperature may cause condensation to develop inside the equipment. Do not turn on the equipment until it has reached room temperature.

### **Cooling**

- Allow the proper amount of cooling downtime (for the X-ray tube to cool down) before the acquisition of the next image.
  - Mode of operation: Continuous operation with intermittent loading - Needs waiting time (at least 60 times the exposure time) before the next exposure begins
  - Column operation time: Max. 2 min. On / 18 min. Off (Ratio 1:9)
- If the temperature inside the tube head reaches 60 °C (140 °F), X-ray exposure will cease, and an error message will be displayed. Normal X-ray capabilities will resume after the generator reaches 58 °C (136.4 °F).
- If the fan (optional) is installed, it automatically operates when the temperature surrounding the tube head reaches the pre-defined level: 40 °C (104 °F). The setpoint temperature is configurable.

### **Turning the equipment on / Adjusting the height of the equipment**

- Do not position the patient near the equipment while it is initiating as the patient could be injured if the equipment malfunctions.
- Ensure that the patient is kept clear of the equipment while adjusting its height.

#### Emergency stop

- If a problem occurs during image acquisition, press the red **Emergency Stop Switch** to immediately stop all moving parts and cut off all power to the equipment. (**Emergency Stop Switch** is located under the bottom of the Vertical Frame. Turn the switch in the direction of the arrow to reboot the equipment.)

#### Trouble-free operation

- Never use this equipment in an environment that is susceptible to explosion.
- Always operate the equipment within a temperature range of 10 °C to 35 °C (50 °F to 95 °F) for the safe operation. Image quality may deteriorate if the equipment is operated outside of this range.
- Always allow the equipment enough time to warm up (while switched on) if it has been exposed to temperatures below 10 °C (50 °F).
- Only perform X-rays of patients if the system is in full working order.
- Always ensure that equipment movement is not obstructed by the patient's clothing, a medical device (such as a wheelchair), or the patient.
- Do not leave the patient unattended around the equipment.
- Remove all radio-controlled devices, mobile phones, etc. from the X-ray room before image acquisition as these objects may cause the equipment to malfunction.

#### Modifying the equipment

- Modifying the equipment in any way which may affect the safety of the operator, patients or other persons is prohibited by law.
- No part of this equipment is serviceable by the operator. A **VATECH** qualified service technician must perform all maintenance and repair of this equipment.
- This product may only be operated with original **VATECH** accessories or third-party accessories expressly approved by **VATECH**.

## 3.2 Electricity-related Safety Precautions



To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.

- Check the status of the power source, PC, and cables before operating the equipment.
- Ensure that **Main Power Switch** is set to off when the equipment is not in use.
- Always disconnect the power supply before cleaning the equipment.
- Always keep electrical cords away from hot appliances or radiators.
- DO NOT place the PC or peripheral equipment connected to the PC near the patient.
- The equipment and PC should be connected to a common protective earth.
- Never overload the equipment's circuit by sharing it with too many appliances.
- Use the same power circuit for the PC and the equipment.

### Combining this equipment with other devices

- Do not connect this equipment to devices that are not designated as a part of the system.
- Do not connect this equipment to a Multiple Portable Socket-Outlet (MPSO) or extension cord which is not provided with the equipment.

### Electromagnetic compatibility

- This X-ray equipment complies with IEC standard 60601-1-2.
- Medical electrical equipment is subject to special Electromagnetic Compatibility (EMC) preventive measures. It must be installed and operated as specified in EMC information.
- If high-voltage systems, radio link systems or MRI systems are located within 5 m of the unit, please observe the specifications stated in the installation requirements.
- Portable Radio Frequency (RF) communications equipment may interfere with medical electrical equipment. Therefore, the use of mobile wireless phones in medical offices or hospital environments must be prohibited.
- For more details, refer to **15.3 Electromagnetic Compatibility (EMC) Information**.
- Please also observe the Electro-Static Discharge (ESD) protective measures described.

#### **Static Discharge**


- Connector pins or sockets bearing ESD warning labels must not be touched or interconnected without observing ESD protective measures.



Electrostatic discharge (ESD)

#### **ESD protective measures include**

- Procedures for preventing electrostatic charge build-up (e.g., temperature control, humidification, conductive floor coverings and non-synthetic clothing)
- Electrostatic discharge of your own body with the frame of the equipment, the protective ground wire or large metallic objects
- Use of the wristband for grounding



### 3.3 Radiation Safety



Since rules and regulations concerning radiation safety differ between countries, it is the responsibility of the owner and operator of this equipment to comply with all applicable rules and regulations concerning radiation safety and protection in his/her area.

- This equipment must be housed inside an X-ray shielded room.
- The operator must remain outside a shielded room during X-ray exposure to protect himself/herself from radiation.
- During imaging; the operator must maintain vocal/visual contact with the patient from outside the shielded area.
- The operator should continuously check the status of the patient and the equipment during imaging.
- The operator should be at least 2 m (6 feet) away from the equipment during imaging.
- The operator must immediately stop imaging if the equipment malfunctions.
- The patient must wear a lead apron with neck and thyroid protection during X-ray exposure.
- Children and pregnant women must consult with a doctor before X-ray exposure.



As a manufacturer of radiology equipment that conforms to stringent protection standards around the world, we guarantee the maximum degree of protection against radiation hazards for our equipment.

### 3.4 Warnings

The following warning statements should be obeyed with the utmost care. Failure to follow these warnings may cause severe damage to the equipment or physical injuries to the patient and the operator.

	<ul style="list-style-type: none"> <li>▪ X-ray equipment is hazardous to the patient and the operator if proper exposure safety measures and operating instructions are not observed.</li> <li>▪ It is essential to read this <b>User Manual</b> carefully and strictly abide by all warnings and cautions stated within it.</li> </ul>
	<ul style="list-style-type: none"> <li>▪ The 3D image should not be used for screening examinations. Each exam must be justified by demonstrating that the benefits outweigh the risk.</li> <li>▪ Where it is likely that evaluation of soft tissues will be required as part of the patient's radiological assessment, conventional medical CT or MR should be used instead of dental cone beam imaging.</li> </ul>
	<ul style="list-style-type: none"> <li>▪ <b>Green16 (PHT-65LHS)</b> system, like other medical equipment, uses high-frequency electrical signals that can interfere with implantable devices such as pacemakers and Implantable Cardioverter Defibrillators (ICDs). If the patient has such an implantable device, you should be aware of any interference in its operation and immediately power off the Dental X-ray system.</li> <li>▪ <b>Green16 (PHT-65LHS)</b> system is designed to withstand the effects of defibrillation. However, when possible, disconnect the Dental X-ray system during defibrillation since a malfunction of the safety controls could otherwise result in electrical burns for the patient.</li> </ul>
	<p>Federal law restricts this device to sale by or on the order of a dentist or with the descriptive designation of any other practitioner licensed by the law of the State in which he practices using or order the use of the device.</p>

### Lasers

- The system incorporates Class 1 laser products. The light localizers used in this product are intended for correct patient positioning and must not be used for any other purpose.
- For maximum safety, advise the patient not to look directly at the laser beam.
- While adjusting the patient, ensure that the laser beam is not directed at the patient's eyes.
- Wavelength: 650 nm, Radiant power: Max. 039 mW



Risk of eye injury!

Do not use this equipment with any other laser sources and do not make any changes to the settings or processes that are described in these operating instructions.

### Cleaning

- Never expose this equipment to liquids, mists or sprays. Exposing this equipment to liquids may cause an electric shock or otherwise damage the system.
- Do not use spray cleaners on this equipment, as this could cause a fire.

### During the Operation

- Never use this equipment in an environment that is susceptible to explosion.
- Do not place flammable materials near this equipment.
- Do not operate the PC while the equipment is operating. Failure to comply with this instruction may result in system malfunction.
- Immediately stop imaging if the equipment malfunctions in any way.
- If a problem occurs during imaging, press the red **Emergency Stop Switch** to immediately stop all moving parts and cut off all power to the equipment's electrical components.
- Never touch the patient while he or she is touching the SIP/SOP connectors.
- The medical electrical equipment or medical electrical system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the medical electrical equipment or medical electrical system should be observed to verify normal operation in the configuration in which it will be used.
- The use of accessories and cables other than those specified, except cables sold by **VATECH** of the medical electrical equipment or medical electrical system as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of EQUIPMENT or SYSTEM.

#### In case of an electrical fire

- Use only fire extinguishers designed for electrical fires to extinguish fires on this equipment.  
  
Liquid extinguishers, such as those which use water, could damage the equipment and cause physical injury.
- Unplug the equipment's power cable before extinguishing any fire.

#### Installation

- To avoid improperly balanced equipment, install the device on a flat surface to maintain stability.
- If the equipment is not stable, property damage and personal injury may occur.
- Do not push or pull the equipment.
- Equipment should only be installed by an authorized technician, complying with proper installation procedures.

#### **NOTICE**

For further details on installation, refer to the **Green16 (PHT-65LHS) Installation Manual**.

#### Security Capabilities

- It is recommended to install and operate **EzDent-i / EasyDent** SW within a secure operating environment that allows only authorized users to access and a system network equipped with Windows built-in firewall, Windows Defender antispyware tools and other commonly used 3<sup>rd</sup> party security tools and application systems.
- The latest updates for anti-virus software and a firewall are recommended.
- The software can be updated by the manufacturer only. Unauthorized software update through a third party, not the manufacturer, is strictly prohibited. For cybersecurity issues related to the software and medical devices, please contact the manufacturer.

## 4. Imaging System Overview

### 4.1 System Components

- **Green16 (PHT-65LHS)** X-ray equipment
- PC system
- Console Software: PANO, CEPH (Optional), CBCT and 3D MODEL Scan
- **EzDent-i / EasyDent:** 2D viewer and patient management software
- **Ez3D-i / Ez3D Plus:** 3D viewer software

### 4.2 Features

- Multi-FOV support: Selectable FOV among 16x9, 12x9, 8x9, 8x5 and 5x5 (cm)
- The multi-imaging solution for Accurate Diagnostics
- Conventional 2D (PANO and CEPH) image acquisition
- 3D scanning for Plaster Cast with FOV 8x9 (cm)
- Touch Screen implemented for easy use
- DICOM (Digital Imaging Communication in Medicine) format supported

### 4.3 Imaging System Options

Configuration	Item	Sensor	
<b>SP</b>	PANO +CBCT	PANO / CBCT	Xmaru1314CF
<b>SC</b>	PANO + CBCT + CEPH	PANO / CBCT	Xmaru1314CF
		CEPH	Xmaru2602CF

## 4.4 Standards and Regulations

### Standards

**Green16 (PHT-65LHS)** is designed and developed to comply with the following international standards and regulations:

- IEC 60601-1, IEC/EN 60601-1-2, IEC 60601-1-3, IEC 60601-1-6, IEC 60601-2-63
- CAN/CSA-C22.2 No. 60601-1:14, CAN/CSA-C22.2 No. 60601-1-3:09, CAN/CSA-C22.2 No. 60601-1-6:11, CAN/CSA-C22.2 No. 60601-2-63:15, CAN/CSA-IEC 62366:15
- ANSI/AAMI ES60601-1:2005 / (R)2012, AND A1:2012, A2:2010 / (R)2012 (Consolidated text - edition 3.1)
- 21 CFR 1020.30, 31, 33
- NEMA Standard publication PS 3.1-3.18, 2008

	This is Class IIb equipment and obtained CE marking in April 2007 for regulations compliance in accordance with the revised European Union's MDD (Medical Devices Directive) 93/42 EEC.
	This equipment received the CSA certification mark by CAN/CSA C22.2 No.601.1 regulations.

### Classifications (IEC 60601-1 6.1)

- The degree of protection against water ingress: Ordinary Equipment: IPX0
- The degree of protection against electric shock: Class 1 equipment, Type B Applied Parts: Temple Supports, Chinrests, and Bites.




## 4.5 Operating Principles

X-ray is emitted when a high voltage is supplied to X-ray tube assembly which frees electrons from the cathode.


They hit anode to produce an X-ray. The machine acquires images by emitting X-ray continuously and rotates on the human tooth at different angles.

Images are acquired, computed and recompiled to reproduce 2D or 3D images.

## 4.6 Imaging System Configuration




### PC Signal Input / Output



No.	Item
1	3D viewer License Key
2	Video output
3	Fiber optic cable (Data in / out)

### 4.7 Equipment Overview




No.	Item	Description
1	X-ray Detector for CEPH (Optional)	Xmaru2602CF for CEPH imaging sensor
2	Nasal Positioner	<ul style="list-style-type: none"> <li>– Positions the patient during CEPH imaging.</li> <li>– The ruler used for reference in an acquired image that is different from the actual size</li> </ul>
3	Ear Rods	Secure the patient’s head during CEPH imaging.
4	Enclosed Component Storage	The place where Bites, Chinrest Assembly, and the other components can be stored.
5	Vertical Frame	Holds the Rotating Unit. Can be controlled with the <b>Column UP/DOWN</b> switch.

No.	Item	Description
6	Rotating Unit	Rotates around the patient's head while the image is being acquired. (Its movement is different according to the scan mode.)
7	X-ray Generator	The vacuum tube where the X-ray is produced.
8	Emergency Stop Switch	Immediately stops the moving parts and cuts off all power to the equipment's electrical components.
9	Main Power Switch	Turns on / off the main power of the equipment.
10	Column UP/DOWN Switch (optional)	Adjusts the height of the Vertical Frame.
11	Stationary Column	Supports the whole part of the equipment.
12	Base (Optional)	Balances the equipment and maintains its safety.
13	LED Lamp	Displays the status of X-ray exposure. - Green: Standby - Yellow: In operation
14	X-ray Detector for PANO / CBCT	Xmaru1314CF for PANO / CBCT imaging sensor
15	Temple Supports	Supports the patient's head by holding the temples. It is used in PANO and CBCT modes.
16	Chinrest	The place to rest the chin.
17	Control Panel (LCD type)	Operates the Horizontal Beam, opens/closes Temple Supports, adjusts the height of the Vertical Frame and prepares for operation when the <b>READY</b> button is pressed. (For the details, refer to <b>4.6.1 Control Panel</b> .) 
18	D-Sub Connector	The input signal port for <b>Column UP/DOWN</b> Switch

### 4.7.1 Control Panel


#### LCD type



No.	Buttons	Description
1	Touch Screen (LCD)	Configures the parameter settings in each imaging mode. For more information, refer to the <b>4.6.2 Touch Screen</b> .
2	Temple Supports OPEN/CLOSE button	Adjusts the Temple Supports for patient positioning.
3	READY / RETURN button	Indicates that imaging is ready after parameter settings and the patient positioning is complete. Initializes the positioning of the Rotating Unit.
4	Horizontal Beam UP/DOWN button	Aligns the Horizontal Beam in PANO mode.
5	Column UP/DOWN button	Moves the Vertical Frame up or down. (For adjusting the height of the Chinrest)

## 4.7.2 Touch Screen

Set options for imaging of each mode by using the Touch Screen. It provides the same function as the PC's Console Software. Touch Screen and Console Software (5.3. Console Software) are interlocked mutually, therefore, indicate the same environment is setting values always.




### NOTICE

The **CEPH** button exists only when the CEPH imaging program is included in the equipment.


### CAUTION

- Do not allow the patient to control Touch Screen. Doing so may cause physical injury to the patient or damage to the equipment.
- Always operate the Touch Screen by pressing it gently with your fingertip.
- Do not use pointed objects such as ballpoint pens or pencils. Doing so may cause damage to the screen.

**PANO Main Screen**




No.	Function	Description
1	Examination mode selection panel	Displays available PANO Examination programs. (In PANO – Orthogonal mode, press <b>UP/DOWN</b> button to scroll through next/previous ROI option.)
2	<b>CONFIRM</b> button	Confirms the settings and moves to the next step.
3	<b>BACK</b> button	Moves back to the modality (PANO / CEPH (Optional) / CBCT / MODEL) selection screen.
4	<b>Settings</b> button	Adjusts Tube voltage, Tube current, Gender / Age group, X-ray Intensity, Examination type, Imaging type, Arch selection, and Image option.
5	Imaging parameter settings information	Displays the current setting information. (Modality, Tube voltage, Tube current, Arch type, Image option, and Pano option.)

**PANO Settings Screen**

No.	Function	Description
1	Examination Program	Selects between Normal and Special.
2	Pano Option	Selects between Normal, Magic PAN (Optional), and Insight PAN (Optional).
3	Arch	Selects available patient's Arch types.
4	Image Option	The default is "High Resolution." "Green" is optional. - When "Green" is enabled, Image Option is selectable between "High Resolution" and "Green." When "Green" is disabled, the Image Options section is invisible. (Image quality: High Resolution > Green)
5	kVp / mA control button	Adjusts Tube voltage (kVp) and Tube current (mA).
6	Patient's Gender / Age group	Selects patient's Gender / Age group.
7	X-ray intensity	Selects X-ray intensity. <b>NOTICE</b> Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft: Soft ≤ Normal ≤ Hard
8	<b>EXIT</b> button	Closes the Settings Screen and moves back to <b>PANO</b> Main Screen.

**CEPH Main Screen**




No.	Function	Description
1	Examination selection panel	Displays available CEPH Examination programs.
2	<b>CONFIRM</b> button	Confirms the settings and moves to the next step
3	<b>BACK</b> button	Moves back to the modality (PANO / CEPH (Optional) / CBCT / MODEL) selection screen.
4	<b>Settings</b> button	Adjusts Tube voltage, Tube current, Gender / Age Group, X-ray Intensity, and Image Option.
5	Imaging parameter settings information	Displays the current setting information. (Modality, Tube voltage, Tube current, and Image option)

**CEPH Settings Screen**




No.	Function	Description
1	kVp / mA control button	Adjusts Tube voltage (kVp) and Tube current (mA).
2	Patient's Gender / Age group	Selects patient's Gender / Age group.
3	X-ray intensity	Selects X-ray intensity.  <div style="border: 1px solid black; padding: 5px; display: inline-block;"><b>NOTICE</b></div> Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal or Soft: Soft ≤ Normal ≤ Hard
4	Image Option	Selects between "High Resolution" and "Green." (Image quality: High Resolution > Green)
5	<b>EXIT</b> button	Closes the Settings Screen and moves back to CEPH Main Screen.

**CBCT Main Screen**




No.	Function	Description
1	FOV selection panel	Displays available FOV modes.
2	<b>CONFIRM</b> button	Confirms the settings and moves to the next step.
3	<b>BACK</b> button	Moves back to the modality (PANO / CEPH (Optional) / CBCT / MODEL) selection screen.
4	<b>Settings</b> button	Adjusts Tube voltage, Tube current, Gender / Age group, X-ray intensity, Image Option, and Voxel Size.
5	Imaging parameter settings information	Displays the current setting information. (Modality, Tube voltage, Tube current, Image option, and Voxel size)

**CBCT Settings Screen**



No.	Function	Description
1	Process Option	Selects between “HD” and “Fast.”
2	Image Option	Selects between “High Resolution” and “Green.” (Image quality: High Resolution > Green)
3	Voxel size	Selects between Standard and Application.
4	kVp / mA control button	Adjusts Tube voltage (kVp) and Tube current (mA).
5	Patient's Gender / Age group	Selects patient's Gender / Age group.
6	X-ray intensity	Selects X-ray intensity.  <div style="border: 1px solid black; padding: 5px; display: inline-block;"><b>NOTICE</b></div> Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft: Soft ≤ Normal ≤ Hard
7	<b>EXIT</b> button	Closes the Settings Screen and moves back to CBCT Main Screen.

**3D MODEL Scan Main Screen**



No.	Function	Description
1	Examination selection panel	Displays available 3D MODEL Scan Examination programs.
2	<b>CONFIRM</b> button	Confirms the settings and moves to the next step
3	<b>BACK</b> button	Moves back to the modality (PANO / CEPH (Optional) / CBCT / MODEL) selection screen.
4	<b>Settings</b> button	Adjusts Tube voltage, Tube current, Gender / Age group, and X-ray intensity.
5	Imaging parameter settings information	Displays the current setting information. (Modality, Tube voltage, and Tube current)

**3D MODEL Scan Settings Screen**

No.	Function	Description
1	kVp / mA control button	Adjusts Tube voltage (kVp) and Tube current (mA).
2	Patient's gender / age group	Selects patient's gender / age group.
3	X-ray intensity	Selects X-ray intensity.
4	<b>EXIT</b> button	Closes the Settings Screen and moves back to 3D MODEL Scan Main Screen.

### 4.7.3 Emergency Stop Switch

During operation, the following emergency situations may occur:

- X-ray emission even after the **Exposure Switch** has been released
- Physical injury to the patient or damage to the equipment
- Other emergency situations

If a problem occurs during image acquisition, press the red **Emergency Stop Switch** to immediately stop the moving parts and cut off all power to the equipment's electrical components. To restart the equipment, turn the **Emergency Stop Switch** clockwise until it pops up.

The **Emergency Stop Switch** is located under the bottom of the Vertical Frame.




### 4.7.4 Exposure Switch

The **Exposure Switch** allows the operator to control image acquisition from outside of the X-ray room.

Press and hold the **Exposure Switch** down until acquisition is completed. Premature release of the **Exposure Switch** will abort image acquisition.

Pressing the **Exposure Switch** activates the LED indicator to turn yellow. This color indicates that the X-ray is being emitted.



#### IMPORTANT

The **Exposure Switch** is detachable. Ensure the **Exposure Switch** cable is not detached from the unit accidentally during the operation.


Keep vocal/visual contact with the patient during exposure. If any problem occurs during exposure, release the **Exposure Switch** immediately.

### 4.7.5 Enclosed Components

The enclosed components can be disassembled and cleaned. All enclosed components that are used to support the patient (the Bite, the Chinrest, the Ear Rods, and the Temple Supports) should be cleaned with ethanol and wiped with clean towels.

Components	Name and Function	Materials
	Normal Bite - For PANO and CBCT normal patients	PC (Polycarbonate)
	Deep Bite Block*	PC (Polycarbonate)
	Special Bite A - For PANO TMJ and Sinus modes - For CBCT TMJ patients	PC (Polycarbonate)
	Special Bite B - For PANO edentulous patients - For CBCT Double scan (Mandible) mode.	PC (Polycarbonate)
	Double Scan Support - For CBCT Double scan (Maxilla) mode.	PC (Polycarbonate)
	Normal Chinrest - For Normal Bite - For CBCT Double scan (Mandible) mode.	ABS (Acrylonitrile butadiene styrene) copolymer
	Special Chinrest - For Special Bite A and Special Bite B	ABS (Acrylonitrile butadiene styrene) copolymer
	Temple Supports (1 set)	PC (Polycarbonate)
	Ear Rods (1 set)	Silicone
	Nasal Positioner Cover - For CEPH	Silicone

## 4. Imaging System Overview

Components	Name and Function	Materials
	Carpus Plate	PC (Polycarbonate)
	Sanitary Vinyl Covers (disposable) for the Bite Block	LDPE (Low-density polyethylene)
	Sanitary Vinyl Covers (disposable) for the Double Scan Support	LDPE (Low-density polyethylene)
	Protractor (1 set) - For positioning the patient's body in CEPH mode.	PC (Polycarbonate)
	Model Scan Jig	ABS (Acrylonitrile butadiene styrene) copolymer

## 5. Imaging Software Overview

Three programs are included in this equipment to acquire, process, and view the image:

- **EzDent-i / EasyDent:** 2D viewer and patient management software
- **Ez3D-i / Ez3D Plus:** 3D viewer software
- Console software: PANO, CEPH (Optional), CBCT and 3D MODEL Scan

### 5.1 PC Specifications (Recommended)

#### NOTICE

- The PC system plays a vital role in image processing and verification. Configure the PC environment to meet the following specifications. If the PC specifications are not met, the image quality can be lower.
- Do not place patients near the equipment and PC.

Item	Specifications
CPU	Intel Xeon E5-1620v3 3.5GHz 2133 4C
RAM	16GB DDR4-2133 Registered RAM or larger
HDD	1TB SATA
Graphics board	NVIDIA GeForce GTX1060 D5 6GB or greater
Ethernet Interface	Integrated Intel I218LM PCIe GbE
Serial Port (RS232)	HP Serial Port Adapter kit
Power Supply	≥ 700 Watts (90% efficient)
Slots	2 PCI Express Gen3 x16 slot 1 PCI Express Gen3 x 8 Slot 1 PCI Express Gen2 x 4 Slot 1 PCI Express Gen2 x 1 Slot 1 PCI Slot
CD/DVD drive	DVD-ROM, DVD+/-RW, Blu-Ray
Monitor	19" 1280x1024 screen resolution
Operating System	Windows 10
Recommended System	HP Z4

## 5.2 EzDent-i / EasyDent

**EzDent-i / EasyDent** is imaging software from **VATECH Co., Ltd.** that manages patient images to make faster and more accurate diagnoses. **EzDent-i / EasyDent**, linked with the console software and 3D viewer, makes it convenient for the operator to use and process necessary images. Various functions enable the acquired images to be processed quickly and conveniently from the console software.

**NOTICE**

Please refer to **5.2.1 Creating a New Patient Record** and **5.2.2 Retrieving Patient Records** and **EzDent-i / EasyDent User Manual** for more information.

**NOTICE****Security Capabilities**

- It is recommended to install and operate **EzDent-i / EasyDent** SW within a secure operating environment that allows only authorized users to access and a system network equipped with Windows built-in firewall, Windows Defender antispysware tools and other commonly used 3<sup>rd</sup> party security tools and application systems.
- The latest updates for anti-virus software and a firewall are recommended.
- The software can be updated by the manufacturer only. Unauthorized software update through a third party, not the manufacturer, is strictly prohibited. For cybersecurity issues related to the software and medical devices, please contact the manufacturer.


### 5.3 Console Software

Use the Console Software to configure the imaging environment according to the mode.

**NOTICE**

- You can set the imaging parameters on either the Touch Screen or the Console Software running on the PC. They are synchronized and display the same environmental settings.)
- To improve program functions, the Console Software may change without notice

The Main Screen of the Console Software consists as follows. Each imaging mode will be described later.



No.	Item	Description
1	Imaging Mode Display	Displays the current imaging mode.
2	Scanning Status and Image Preview window	Shows image acquisition progression and display a preview of the images acquired.
3	Patient Information	Displays information on the selected patient.

No.	Item	Description																				
	window																					
4	Guide Message window	Displays various text instructions for the operator.																				
5	DAP, Scan Time and Exposure Time Display window	Displays estimated DAP (Dose Area Product), scan time and exposure time after exposure parameter settings are completed.																				
6	Tube Voltage and Current Adjustment	<p>If the patient is selected in <b>EzDent-i / EasyDent</b>, the default kVp / mA according to the patient's information (gender/age) is displayed. This tool adjusts the kVp and mA values and controls the power of the X-ray to improve image quality. If necessary, adjust the kVp and mA values manually using the arrows.</p> <div style="border: 1px dashed black; padding: 5px;"> <div style="display: flex; align-items: center;"> <div style="background-color: #0056b3; color: white; padding: 5px; font-weight: bold; margin-right: 10px;">NOTICE</div> <div> <p>For the tube voltage and its correspondence with the selected patient, refer to <b>16.1 Recommended X-ray Exposure Table</b>.</p> </div> </div> </div>																				
7	Patient's gender/age group	<p>Displays the current patient's gender/age group as entered in <b>EzDent-i / EasyDent</b>'s patient information fields. If necessary, the gender/ age group can be manually selected.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Gender / Age Group</th> <th>VATECH's Standard</th> </tr> </thead> <tbody> <tr> <td colspan="2" style="text-align: center;">Child</td> <td style="text-align: center;">2 ~ 12 years of age</td> </tr> <tr> <td rowspan="2" style="text-align: center;">Adult</td> <td style="text-align: center;">Man</td> <td rowspan="2" style="text-align: center;">&gt; 12 years of age</td> </tr> <tr> <td style="text-align: center;">Woman</td> </tr> </tbody> </table>	Gender / Age Group		VATECH's Standard	Child		2 ~ 12 years of age	Adult	Man	> 12 years of age	Woman										
Gender / Age Group		VATECH's Standard																				
Child		2 ~ 12 years of age																				
Adult	Man	> 12 years of age																				
	Woman																					
8	X-ray intensity	<p>Selects X-ray intensity.</p> <div style="border: 1px dashed black; padding: 5px;"> <div style="display: flex; align-items: center;"> <div style="background-color: #0056b3; color: white; padding: 5px; font-weight: bold; margin-right: 10px;">NOTICE</div> <div> <p>Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft: Soft ≤ Normal ≤ Hard</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Age Group</th> <th>Average Head Circumference (cm)</th> <th>Range (cm)</th> <th>X-ray Intensity</th> </tr> </thead> <tbody> <tr> <td rowspan="3" style="text-align: center;">Child</td> <td rowspan="3" style="text-align: center;">53±3</td> <td style="text-align: center;">&gt;53±3</td> <td style="text-align: center;">Hard</td> </tr> <tr> <td style="text-align: center;">53±3</td> <td style="text-align: center;">Normal</td> </tr> <tr> <td style="text-align: center;">&lt;53±3</td> <td style="text-align: center;">Soft</td> </tr> <tr> <td rowspan="3" style="text-align: center;">Adult</td> <td rowspan="3" style="text-align: center;">56±3</td> <td style="text-align: center;">&gt;56±3</td> <td style="text-align: center;">Hard</td> </tr> <tr> <td style="text-align: center;">56±3</td> <td style="text-align: center;">Normal</td> </tr> <tr> <td style="text-align: center;">&lt;56±3</td> <td style="text-align: center;">Soft</td> </tr> </tbody> </table> </div> </div> </div>	Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity	Child	53±3	>53±3	Hard	53±3	Normal	<53±3	Soft	Adult	56±3	>56±3	Hard	56±3	Normal	<56±3	Soft
Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity																			
Child	53±3	>53±3	Hard																			
		53±3	Normal																			
		<53±3	Soft																			
Adult	56±3	>56±3	Hard																			
		56±3	Normal																			
		<56±3	Soft																			
9	<b>CONFIRM / READY</b> button	<p><b>CONFIRM</b></p> <p>Applies the selected settings and moves to the next step. (Exposure parameter setting and patient positioning &gt; Ready for exposure)</p>																				

No.	Item	Description
		<div data-bbox="605 336 817 392" style="border: 1px solid black; background-color: #0056b3; color: white; padding: 5px; text-align: center; font-weight: bold;">NOTICE</div> <p data-bbox="838 282 1208 446">When you click the <b>CONFIRM</b> button, estimated DAP (Dose Area Product), scan time and exposure time would be displayed DAP, Scan Time and Exposure Time Display window.</p> <div data-bbox="587 459 714 498" style="background-color: #cccccc; padding: 2px;"><b>READY</b></div> <p data-bbox="587 508 1208 614">It is activated when you click the <b>CONFIRM</b> button after the patient positioning is completed. Click the button when all aspects of preparation are completed for image acquisition.</p>
10	Imaging parameters configuration panel	Selects the imaging parameters for each mode: PANO, CEPH (Optional), CBCT and 3D MODEL Scan
11	<b>Modality Selection</b> button	Returns to Modality Selection (PANO / CEPH (Optional) / CBCT / MODEL) screen.
12	Test Rotation button	<p data-bbox="587 832 1192 977">Switches to the Test mode to check if any part of the patient's body is reached to the surface of the equipment before the actual exposure. This function is applied to PANO and CBCT modalities only.</p> <p data-bbox="587 987 883 1012">To change to the test mode,</p> <ol data-bbox="642 1022 1208 1234" style="list-style-type: none"> <li>1. Align the patient to the equipment. (For the details, refer to the "Positioning the Patient" section of each modality chapter.)</li> <li>2. Select a modality.</li> <li>3. Click the CONFIRM button.</li> <li>4. Click the Test Rotation button. Then, "ON" on the Test Rotation button is changed to "OFF."</li> </ol> <div data-bbox="707 1244 1044 1344" style="text-align: center;">  </div> <ul data-bbox="696 1354 1167 1466" style="list-style-type: none"> <li>● To initiate test rotation, press the <b>BEAM ON/OFF</b> button on the Control Panel.</li> <li>● To finish the test mode, click the <b>Test Rotation</b> button or <b>READY</b> button.</li> </ul> <div data-bbox="605 1495 817 1551" style="border: 1px solid black; background-color: #0056b3; color: white; padding: 5px; text-align: center; font-weight: bold;">NOTICE</div> <p data-bbox="838 1495 1188 1551">This function is applied to PANO and CBCT modalities only.</p>

No.	Item	Description
13	Phantom button	<p>This function is used to acquire Phantom images.</p> <hr/> <div style="display: flex; align-items: center;"> <div style="border: 1px solid black; background-color: #0056b3; color: white; padding: 5px; margin-right: 10px;"><b>NOTICE</b></div> <div> <p><b>Image acquisition using the Phantom Jig:</b></p> <ol style="list-style-type: none"> <li>1. Click the <b>Phantom</b> button.</li> <li>2. Select the Modality and click the <b>Capture</b> button.</li> <li>3. Check the parameters displayed in the main GUI window and align the Phantom Jig, and then click the <b>READY</b> button.</li> <li>4. Press and hold down the <b>Exposure Switch</b>.</li> </ol> </div> </div>
14	speaker volume button	<p>This button is used to adjust the speaker volume. Clicking on the speaker icon brings up the volume control bar, and you can adjust the volume by clicking and moving the volume control bar with your mouse. After moving the bar, release the mouse to play the current volume and save the current volume.</p>
15	Manual Reconstruction button	<p>Reconstructs the image manually when automatic image reconstruction fails:            Select a Modality after clicking this button. &gt; Click the <b>Search</b> button. &gt; Select an image to reconstruct. &gt; Click the <b>Reconstruction</b> button.</p>
16	Laser Beam ON/OFF button	<p>Turns the Laser Beam on or off for patient positioning. Enabled when the <b>CONFIRM</b> button is clicked after the imaging conditions are configured.</p>
17	Settings button	<p>Displays and sets various equipment-related parameters, including language, automatic save, DAP display unit, etc.</p>
18	EXIT button	<p>Exits the console software.</p>
19	X-ray indicator	<p>The radiation mark turns yellow and “X-RAY” changes to “X-RAY ON.”</p> <div style="display: flex; align-items: center;">  <div style="background-color: #f0f0f0; padding: 5px; border: 1px solid #ccc;">X-RAY ON</div> </div>
20	Version Information	<p>Displays the Console Software version.</p>

## 6. Getting Started

### 6.1 Turning on the Equipment



- Do not place the patient close to the equipment when it's being turned on. Doing so may cause physical injury to the patient and damage to the equipment.
- Do not operate the PC while the equipment is in operation. Doing so may cause an error in the equipment.



- The extreme fluctuation of temperature may cause condensation inside the equipment. Do not switch on the equipment until it has reached normal room temperature.
- Rebooting the equipment: After turning it off, wait for approx. — 20 seconds before turning it on again.
- Warm-up the equipment for at least 5 minutes before the operation. For the best image quality, it is recommended to have a warm-up phase for more than 30 minutes.

#### IMPORTANT

If the equipment has not been used for a long time, please let it have enough time to be warmed up. It extends the life of the X-ray tube.

The imaging system mainly consists of the imaging equipment and the PC.

Before turning on the equipment, please confirm that the equipment and PC have been installed correctly.

- Turn on the PC.
- Press the **Main Power Switch** that is located under the bottom of the Vertical Frame to turn on the equipment.



#### NOTICE

**Main Power Switch** isolates its circuits electrically from the supply mains on all poles simultaneously.

- Make sure that the green LED light at the top of the equipment is on.

## 6.2 Running the Image Viewer (EzDent-i / EasyDent)

The Imaging Program is interfaced with **EzDent-i / EasyDent**, and the user can analyze the image acquired from the Console Software easily and rapidly. On your desktop, double-click the **EzDent-i / EasyDent** icon. The **EzDent-i / EasyDent** main window will be displayed.

### NOTICE

For further details on this subject, refer to the **EzDent-i / EasyDent User Manual**.

### NOTICE

#### Security Capabilities

- It is recommended to install and operate **EzDent-i / EasyDent** SW within a secure operating environment that allows only authorized users to access and a system network equipped with Windows built-in firewall, Windows Defender antispysware tools and other commonly used 3<sup>rd</sup> party security tools and application systems.
- The latest updates for anti-virus software and a firewall are recommended.
- The software can be updated by the manufacturer only. Unauthorized software update through a third party, not the manufacturer, is strictly prohibited. For cybersecurity issues related to the software and medical devices, please contact the manufacturer.

### NOTICE

For **Green16 (PHT-65LHS)** dental computed tomography X-ray system, both 3D viewer (**Ez3D-i / Ez3D Plus**) and Console Software are being accessed through 2D viewer (**EzDent-i / EasyDent**) SW. 3D viewer and Console Software do not have an image storage capacity of their own, and both programs will not be able to keep patient information.

## 6.2.1 Creating a New Patient Record

To create a new patient record, follow the procedure outlined below:


### EzDent-i

1. Click the **PATIENT** tab and click the **Add Patient** icon from the main GUI window.

2. Enter the required patient information. **Chart Number**, **E-Mail address**, **First Name**, and **Last Name** are required fields that must be filled in. (The Chart Number is filled in automatically.)
3. Click the **Add** button to save the patient record.

## EasyDent

1. Click the **Patient** icon in the upper left corner of the **EasyDent's** main GUI window.



2. Enter the required patient information. **Chart Number**, **First Name**, and **Last Name** are required fields that must be filled in. All other fields are optional, but it is recommended that they are filled in.

3. Click the **Add** button to save the patient record.

## 6.2.2 Retrieving Patient Records

You can search through the patient database using a patient's Chart Number, First Name, or Last Name.

### EzDent-i

1. Enter the Name or Chart Number of the patient to be searched on the **Patient Search** panel and then click the **Search** button. The information on the patient that fits the search condition appears.




Double-click the Keyboard icon to display the virtual keyboard. You may search for patient information using the virtual keyboard.

**NOTICE**



2. Double-click the patient information to see more details about the patient as shown below.




## EasyDent

1. On the **Patient Information** pane, double-click the Chart No., First name, or Last name of the patient and the virtual keyboard will pop up.



2. Enter the **Chart No.**, **First name**, or **Last name** of the patient by clicking the mouse on the virtual keyboard and click **Enter** key. (The physical keyboard can be used to do the same job)
3. Patient information can be displayed on the **Patient Information** panel and **Patient List**.



## 6.3 Initiating the Console Software

### NOTICE


For a new patient, the patient information needs to be registered first.

### EzDent-i

1. Search and select the patient to be captured.



2. Click the **ACQUISITION** tab and select the imaging mode (CT, Panorama, or Cephalo).



3. The Main Screen for the selected mode appears. From the Main Screen, you can configure the imaging parameter settings before acquiring an image.

### NOTICE

Refer to the following **chapters (7 ~ 10)** for more information on image acquisition.

## EasyDent

1. First, click the patient information in the **Patient List**, and click an imaging modality button to select on the upper left corner.



2. The Main Screen for the selected mode appears. From the Main Screen, you can configure the imaging parameter settings before acquiring an image.

### **NOTICE**

Refer to the following **chapters (7 ~ 10)** for more information on image acquisition.

## 7. Acquiring PANO Images

### 7.1 PANO Imaging Program Overview

- **Result Images**





It provides conventional 2D panoramic images.

- **Image Acquisition Method**

It reconstructs U-shaped arch data to a single 2D image utilizing multiple images taken with the X-ray beam scanning specific oral & maxillofacial regions at different angles.

- **Examination Programs**




It is classified as below based on the ROI (Region of Interest).



Examination Type	Arch Selection	ROI	Example
PANO Examination	Narrow Normal Wide Child Orthogonal	Standard	
		Right	
		Front	
		Left	
	Orthogonal	Bitewing*	
		Bitewing Incisor* (Optional)	



Examination Type	Arch Selection	ROI	Example
		Bitewing Right*	
		Bitewing Left*	
SPECIAL Examination	N/A	TMJ LAT Open	
		TMJ LAT Close	
		TMJ PA Open (Optional)	
		TMJ PA Close (Optional)	
		Sinus LAT (Optional)	
		Sinus PA	

\* Bitewing imaging mode is activated only when Orthogonal is selected in Arch Selection.

■ Main Imaging Programs

Examination Type	Arch Selection	ROI	Description & Sample Image
PANO Examination	Narrow	Standard	<p>A panoramic imaging mode for patients with a V-shaped arch trajectory. (Typically for some females)</p> 
	Normal	Standard	<p>A panoramic imaging mode for adult patients with a typical arch trajectory.</p> 
	Wide	Standard	<p>A panoramic imaging mode for the patients with a square-shaped arch trajectory. (Typically for some males)</p> 
	Child	Standard	<p>A panoramic imaging mode for child trajectory. (Less X-ray exposure than the Normal mode by approximately more than 40%)</p>

Examination Type	Arch Selection	ROI	Description & Sample Image
			
	Orthogonal	Standard	<p>A panoramic imaging mode to minimize the overlapped region of the teeth from the X-ray exposure which is beamed perpendicularly between teeth.</p> 
		Bitewing** (Bitewing Incisor mode is Optional)	<p>A panoramic imaging mode to acquire an image only for the region of interest through the orthogonal trajectory. (Pros: less X-ray exposure than the Normal mode. / Cons: TMJ and some parts of the maxillary sinus cannot be acquired.)</p> 
SPECIAL Examination	N/A	TMJ LAT Open / Close	An imaging mode to acquire a lateral image of the TMJ, in which the X-ray

Examination Type	Arch Selection	ROI	Description & Sample Image
			<p>beam is directed on the lateral TMJ region. (TMJ Open and Close)</p> 
		<p>TMJ PA Open / Close (Optional)</p>	<p>An imaging mode to acquire a TMJ image, in which the X-ray beam is directed on the frontal TMJ, with the patient's mouth open fully and close.</p> 
		<p>Sinus LAT (Optional)</p>	<p>A special imaging mode to acquire a Sinus image, in which an X-ray beam is directed on the lateral region of the maxillary sinus.</p> 
		<p>Sinus PA</p>	<p>A special imaging mode to acquire a Sinus image, in which an X-ray beam is directed on the frontal region of the maxillary sinus.</p> 


\*\* Bitewing imaging mode is activated only when Orthogonal is selected in Arch Selection.

## 7.2 Configuring Exposure Parameters


To acquire PANO Images, **6. Getting Started** must be completed first. If not, you must return to the **6. Getting Started** and finish the step first.

### NOTICE

You can set the imaging parameters on either the Touch Screen or the Console Software running on the PC. They are synchronized and display the same environmental settings.



1. Click the **PANO** button on the Main Screen.



**NOTICE** The **CEPH** button exists only when the CEPH imaging program is included in the equipment.

2. Select a Pano Option. (On Touch Screen, click **Settings** button before selecting options.)




Mode	Description
Normal	- Provides a normal panoramic image.
Magic PAN (Optional)	- Provides a single optimal panoramic image having multiple focal images combined. - Minimizes the difference in the quality of images varied according to the patient's positioning and arch shape.
Insight PAN (Optional)	- Provides multiple panoramic images having different focal planes along with a typical panoramic image together. - Enables detailed verification of images in-depth direction.

**NOTICE** When "Insight PAN" is selected, Image Options are disabled.

3. Select an Image Option.



<Console Software>



<Touch Screen>

Mode	Description
High Resolution	High-Resolution image
Green	Normal quality image


**NOTICE**

The default is "High Resolution." "Green" is optional. When "Green" is enabled, Image Option is selectable between "High Resolution" and "Green." When "Green" is disabled, the Image Options section is invisible.

4. Make an Arch Selection.



<Console Software>



<Touch Screen>

Arch Selection	Description
Narrow	Panoramic image of V-shaped palatal arches (a small number of adult females)
Normal	Panoramic image of normal adult palatal arches
Wide	Panoramic image of square-shaped palatal arches (some number of adult males)
Child	Panoramic image of child palatal arches, approximately more than 40% less X-ray dose than in Normal mode.
Orthogonal	<p>A panoramic image where the x-ray angle enters vertically in between the teeth, so overlapping images is minimized.</p> <p><b>NOTICE</b></p> <p>If Orthogonal Arch is selected, Bitewing examinations (Bitewing, Bitewing Incisor (Optional), Bitewing Right, Bitewing Left) are activated.</p>

- Select an Examination Program in the Pano Examination or Special Examination panel.



<Console Software>



<Touch Screen>

**NOTICE**

- To activate Bitewing examination options- Bitewing, Bitewing Incisor (Optional), Bitewing Right, Bitewing Left, select Orthogonal Arch in the Arch Selection panel.



- When a Special Examination option is clicked, the "PANO Examination" panel is disabled. If you want to select a PANO Examination option, please conduct Arch selection again.
- For more information about the Examination Program, refer to the **7.1 PANO Imaging Program Overview**.

- The Gender / Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.



<Console Software>



<Touch Screen>

**NOTICE**

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

7. Select X-ray intensity.



<Console Software>



<Touch Screen>

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft:

Soft ≤ Normal ≤ Hard


**NOTICE**

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

8. The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ±1 kVp and ±1 mA respectively.



<Console Software>



<Touch Screen>

9. Click the **CONFIRM** button when the exposure parameter setting is completed.



<Console Software>



<Touch Screen>

**NOTICE**

When you click **CONFIRM** button,

- The Rotating Unit will move to its initial scanning position.
- The Vertical Beam will be activated to make patient positioning easier.
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below the Patient Information window.

DAP  
127.334307 mGy x cm<sup>2</sup>

Scan-time  
13.5 Sec

Exposure-time  
13.5 Sec

10. Guide the patient to the equipment.

## 7.3 Patient Positioning




- Have patience (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation.
- Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

### IMPORTANT

- Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition.
- Metal implants or bridges may reduce the quality of the images.
- Be sure to adjust the laser beam correctly. Otherwise, the quality of images can be lower due to ghost images or expansion/reduction of the images.

### Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hairpins, braces, false teeth, etc.). Metal objects may induce ghost images and lower image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button or switch option to adjust the equipment to match the height of the patient.




< Control Panel – LCD type >

### 7.3.1 PANO Examination Mode (Standard / Right / Left / Front / Orthogonal)

#### Normal Patient Positioning

1. Insert the Normal Bite into the Normal Chinrest and cover it with a Sanitary Vinyl Cover.




- The Sanitary Vinyl Cover is for single use only. It should be replaced for each patient. Be sure to use the approved vinyl cover.
- Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.

2. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.




3. Guide the patient to the inside of the equipment.
4. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.
5. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.

6. Let the patient bite the Bite along its grooves with his/her front teeth.




7. Let the patient maintain the posture as follows:

- Close the mouth.
- Place the tongue to the roof of the mouth.
- Close the eyes.



### **Edentulous Patient Positioning**

1. Remove the **Normal Chinrest** and insert the **Special Chinrest** into the equipment.
2. Insert the **Special Bite B** into the **Special Chinrest**.



Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.

3. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.



4. Guide the patient to the equipment.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.
6. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.

7. Let the patient maintain the posture as follows:
- Close the mouth.
  - Place the tongue to the roof of the mouth.
  - Close the eyes.



### Laser Beam Aligning



**WARNING**


Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.




**CAUTION**

If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the laser beam properly.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It's to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** button on the control panel to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.



3. Direct the patient to smile and align the Canine Teeth Beam to the center of the canines. Use the Canine Teeth Beam Lever to adjust the position of the beam.



### Finishing Patient Positioning

1. After checking the positions of the patient and the laser beam, click the **Temple Supports OPEN/CLOSE** button on the control panel to prevent the patient's head from moving.



Make sure that the Temple Supports are in the **CLOSE** position before clicking the **READY** button.

2. Click the **READY** button. X-ray exposure has not started yet.
3. Now go to **7.4 X-ray Exposure** to start the exposure.

### 7.3.2 SPECIAL Examination Mode (TMJ / Sinus)

#### <TMJ Open Mode (LAT / PA)>

The TMJ Close image can be acquired after the TMJ Open image is acquired.


#### NOTICE

##### Steps for TMJ Mode

Patient positioning for TMJ Open > Laser Beam Aligning > X-ray Exposure > Patient positioning for TMJ Close > Laser Beam Aligning > X-ray Exposure

#### Patient Positioning

1. Remove the **Normal Chinrest** and insert the **Special Chinrest** into the equipment.
2. Insert the **Special Bite A** into the **Special Chinrest**



#### CAUTION

Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.

3. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.




4. Guide the patient to the equipment.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.

6. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.
7. Guide the patient to press the base of the nose (acanthion point) against the Chinrest and tilt the head forward about 5°. At this point, make sure the patient's jaw does not touch the equipment.

**IMPORTANT**

- If the jaw touches the equipment, it is difficult to maintain the proper position to get good images.
- Be careful; the patient does not touch the equipment with his/her jaw.

8. Let the patient maintain the posture as follows:
  - Open the mouth.
  - Place the tongue to the roof of the mouth.
  - Close the eyes.

**NOTICE**

- As shown in the picture, the support unit of the integrated Chinrest should touch the patient's acanthion point.
- Ask the patient to maintain his/her position until the operation is completed.

### Laser Beam Aligning




Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.



If the laser beam is not correctly positioned there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the laser beam properly.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It's to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** button on the control panel to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.



### Finishing Patient Positioning

1. After checking the positions of the patient and the laser beam, click the **Temple Supports OPEN/CLOSE** button on the control panel to prevent the patient's head from moving.



Make sure that the Temple Supports are in the CLOSE position before clicking the **READY** button.

2. Click the **READY** button. X-ray exposure has not started yet.
3. Now go to **7.4 X-ray Exposure** to start the exposure.

**<TMJ Close Mode (LAT / PA)>**

The TMJ Close image can be acquired after the TMJ Open image is acquired.

**NOTICE****Steps for TMJ Mode**

Patient positioning for TMJ Open > Laser Beam Aligning > X-ray Exposure > Patient positioning for TMJ Close > Laser Beam Aligning > X-ray Exposure

**Patient Positioning**

1. "Do you want to capture a TMJ Close image?" message will appear when the TMJ Open mode is completed. Press/Click **OK** button to begin TMJ Close mode.




<Console Software>



<Touch Screen>

2. Guide the patient to the equipment.
3. Guide the patient to place the base of his/her nose (acanthion point) against the Chinrest and bend the head forward about 5°.
4. Let the patient maintain the posture as follows:
  - Close the mouth.
  - Place the tongue to the roof of the mouth.
  - Close the eyes.



### **NOTICE**

- As shown in the picture, the support unit of the integrated Chinrest should touch the patient's acanthion point.
- Let the patient maintain his/her position until the operation is completed.

### **Laser Beam Aligning**

This is the same as the one for TMJ Open mode.


### **Finishing Patient Positioning**

This is the same as the one for TMJ Open mode.

## &lt;Sinus Mode (LAT / PA)&gt;

**Patient Positioning**

1. Remove the **Normal Chinrest** and insert the **Special Chinrest** into the equipment.
2. Insert the **Special Bite A** into the **Special Chinrest**



Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.

3. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.




4. Guide the patient to the equipment.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.
6. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.
7. Guide the patient to press the base of the nose (acanthion point) against the Chinrest and tilt the head forward about 5°. At this point, make sure the patient's jaw does not touch the equipment.

**IMPORTANT**

- If the jaw touches the equipment, it is difficult to maintain the proper position to get good images.
- Be careful; the patient does not touch the equipment with his/her jaw.

8. Let the patient maintain the posture as follows:

- Close the mouth.
- Place the tongue to the roof of the mouth.
- Close the eyes.



### **NOTICE**

- As shown in the picture, the support unit of the integrated Chinrest should touch the patient's acanthion point.
- Ask the patient to maintain his/her position until the operation is completed.

### Laser Beam Aligning



Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.



If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the Laser Beam properly.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It's to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** button on the control panel to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.

### Finishing Patient Positioning

1. After checking the positions of the patient and the Laser Beam, click the **Temple Supports OPEN/CLOSE** button on the control panel to prevent the patient's head from moving.



Make sure that the Temple Supports are in the CLOSE position before clicking the **READY** button.

2. Click the **READY** button. X-ray exposure has not started yet.
3. Now go to **7.4 X-ray Exposure** to start the exposure.


## 7.4 X-ray Exposure

	<ul style="list-style-type: none"> <li>▪ If an emergency occurs during image acquisition, release the <b>Exposure Switch</b> to cease X-ray emission.</li> <li>▪ The operator shall always observe the X-ray safety regulations applicable to his/her area during the operation of this equipment.</li> </ul>
	<ul style="list-style-type: none"> <li>▪ The operator must always keep vocal/visual contact with the patient during the image acquisition process.</li> <li>▪ Do not operate the PC during exposure. Doing so may cause the system to malfunction.</li> </ul>
	<ul style="list-style-type: none"> <li>▪ Let the patient close the eyes during the operation.</li> <li>▪ To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also, don't let the patient move until the Temple Supports are open.</li> </ul>

1. Get out of the X-ray room and close the door.

	<p>The operator must always keep vocal/visual contact with the patient during image acquisition.</p>
---	--

2. Press and hold down the **Exposure Switch** until image acquisition is completed.



	<p>The image appears on the screen.</p>
	<p>During X-ray exposure, the status appears as follows.</p> <ul style="list-style-type: none"> <li>▪ The LED light of the <b>Exposure Switch</b> turns yellow.</li> <li>▪ The LED light on the top of the equipment turns yellow.</li> <li>▪ An alert sound comes out to indicate that X-ray emission is currently underway.</li> <li>▪ On Console Software, the radiation mark turns yellow and "X-ray" changes to "X-RAY ON."</li> </ul> <div data-bbox="559 1547 779 1605" style="text-align: center;">  </div>

3. Release the **Exposure Switch** when "Image capturing is completed" message appears on the screen.

## 7.5 Finishing the Scan

1. Open the Temple Supports and guide the patient out of the equipment.
2. For Normal Bite, remove the Sanitary Vinyl Cover from the Bite.
3. Press the **READY** button to bring the Rotating Unit back to its initial position.

## 7.6 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i / EasyDent**.

### NOTICE

Refer to the **EzDent-i / EasyDent User Manual** for more information.

1. The images are transferred to **EzDent-i / EasyDent** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.

Left blank intentionally

## 8. Acquiring CEPH Images (Optional)

### 8.1 CEPH Imaging Program Overview

#### ■ Result Images

It provides conventional 2D cephalometric images.

#### ■ Image Acquisition Method

It acquires multiple images by scanning the specific oral & maxillofacial regions with the linear movement of the linear detector and reconstructs them to a single 2D image through computer calculations.

#### ■ Examination Programs

It is classified as below based on the ROI (Region of Interest).

Examination Area	Description	Position
Lateral / Full Lateral (Optional)	<ul style="list-style-type: none"> <li>Used to study craniofacial disease, trauma, and congenital malformation and examine the soft tissue in the otorhinolaryngological area, the sinus, and the hard palate.</li> <li>Measures the angles formed by the connecting lines between the cranial measurement points to further assess the growth of the facial region. It's widely used in Orthodontics and Oral and Maxillofacial Surgery.</li> </ul>	 <p>&lt;Lateral&gt;</p>  <p>&lt;Full Lateral&gt;</p>
PA	<ul style="list-style-type: none"> <li>The radiation is directed from the posterior of the skull to the anterior.</li> <li>It is used to examine cranial diseases, trauma, and congenital malformations.</li> <li>It is used to assess the growth of the lateral side of the face. It is also used to examine the ramus mandibulae, the posterior region of the third-largest molar in the lower jaw, the sidewall of the maxillary sinus, the frontal sinus, antrum ethmoidal, olfactory pits and optic disc pits.</li> <li>Measures the angles formed by the connecting lines between the cranial measurement points to further assess the growth of the facial region. It is widely used in Orthodontics and Oral and Maxillofacial Surgery.</li> </ul>	 <p>&lt;PA&gt;</p>

## 8. Acquiring CEPH Images (Optional)


Examination Area	Description	Position
SMV	<ul style="list-style-type: none"> <li>Used to study the base of the skull, horizontal angulation of the mandibular condylar axis, the sphenoid sinus, the curvature of the lower jaw, the side wall of the maxillary sinus, and zygomatic arch fractures. Also used to study the inner and outer alar plates and holes at the base of the skull.</li> </ul>	 <p data-bbox="1026 517 1104 546">&lt;SMV&gt;</p>
Waters' view	<ul style="list-style-type: none"> <li>Used to study the frontal sinus, the antrum ethmoidal, the optic disc pit, the frontozygomatic suture, the nasal cavity, the coronoid process between the upper jaw and the zygomatic arch.</li> </ul>	 <p data-bbox="989 772 1141 801">&lt;Waters' view&gt;</p>
Carpus	<ul style="list-style-type: none"> <li>Used to assess hand bone age to compare the changes in the skull.</li> </ul>	 <p data-bbox="1016 989 1114 1018">&lt;Carpus&gt;</p>

## 8.2 Configuring Exposure Parameters

To acquire CEPH images, **6. Getting Started** must be completed first.


### NOTICE

You can set the imaging parameters on either the Touch Screen or the Console Software running on the PC. They are synchronized and display the same environmental settings.



## 8. Acquiring CEPH Images (Optional)

1. Click the **CEPH** button on the Main Screen.




<Console Software>

<Touch Screen>


### NOTICE

The **CEPH** button exists only when the CEPH imaging program is included in the equipment.

2. Select an Image Option. (On Touch Screen, click **Settings** button before selecting options.)




<Console Software>



<Touch Screen>

3. Select an examination program in the CEPH Examination panel.



<Console Software>



<Touch Screen>

4. The Gender / Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.



<Console Software>



<Touch Screen>

**NOTICE**

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

5. Select X-ray intensity.



<Console Software>



<Touch Screen>

**NOTICE**

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft:

Soft ≤ Normal ≤ Hard

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

## 8. Acquiring CEPH Images (Optional)

- The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by  $\pm 1$  kVp and  $\pm 1$  mA respectively.



<Console Software>




<Touch Screen>

- Click the **CONFIRM** button when the exposure parameter setting is completed.



<Console Software>



<Touch Screen>

### NOTICE

When you click **CONFIRM** button,

- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below the Patient Information window.

DAP  
127.334307 mGy x cm<sup>2</sup>

Scan-time  
13.5 Sec

Exposure-time  
13.5 Sec

- Guide the patient to the equipment.

### 8.3 Patient Positioning



- Have patience (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation.
- Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.




Ensure that the Nasal Positioner is left unfolded, before adjusting the Ear Rods in the proper direction.

**IMPORTANT**

- Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition.
- Metal implants or bridges may reduce the quality of the images.

#### Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hairpins, braces, false teeth, etc.). Metal objects may induce ghost images and lower image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button or switch option to adjust the equipment to match the height of the patient.



< Control Panel – LCD type >

### 8.3.1 Lateral / Full Lateral (Optional) Mode

#### NOTICE


Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition.

#### Patient Positioning

1. Turn the Nasal Positioner to the **Lateral** mode Positioning Marker as shown below.



2. Leave enough space between the Ear Rods.



#### NOTICE


Use the Porion Position Reference Indicator that appears in the acquired image to quickly confirm the location of Porion.

3. Guide the patient to the CEPH unit.
4. Direct the patient to relax his/her neck and shoulders and stand upright.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.


#### WARNING

After adjusting the height of the column, align the Ear Rods and Nasal Positioner to the patient.

- Align the Ear Rods into the patient's ears properly so that the head does not move during the operation. Moreover, aligning the Nasal Positioner with the patient's nasion by adjusting its height.



- Align horizontally, so the patient's Frankfurt Line is parallel with the floor.



- Direct the patient to swallow first before closing the mouth and to remain in his/her current position until image acquisition is completed.
- Click the **READY** button. The x-ray exposure has not started yet.
- Now go to **8.4 X-ray Exposure** to start the exposure.


### 8.3.2 PA Mode

#### Patient Positioning

1. Turn the Nasal Positioner to the **PA / Waters' view / Carpus** mode Positioning Marker as shown below.



2. Fold the Nasal Positioner up. The Nasal Positioner is not used in PA mode.
3. Guide the patient to the CEPH unit.
4. Ask the patient to stand upright facing the sensor. Make sure that the patient's shoulders are level and that his/her neck is relaxed.




5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.

**WARNING**

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, correctly align the Ear Rods to the patient's ears, so his/her head does not move.
7. Align horizontally, so the patient's Frankfurt Line is parallel with the floor.



8. Direct the patient to swallow first before closing his/her mouth and to remain in his/her current position until image acquisition is completed.
9. Click the **READY** button. The x-ray exposure has not started yet.
10. Now go to **8.4 X-ray Exposure** to start the exposure.


### 8.3.3 SMV Mode

#### Patient Positioning

1. Turn the Nasal Positioner to the **SMV** mode Positioning Marker as shown below.



2. Fold the Nasal Positioner up. The Nasal Positioner is not used in SMV mode.
3. Guide the patient to the CEPH unit.
4. Guide the patient to face the X-ray tube and stand upright.




5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.

**WARNING**

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, correctly align the Ear Rods to the patient's ears, so his/her head does not move.
7. Carefully tilt the patient's head back and adjust so his/her Frankfurt Line is vertical with the floor.
8. Direct the patient to swallow first before closing his/her mouth and to remain in his/her current position until image acquisition is completed.



9. Click the **READY** button. The x-ray exposure has not started yet.
10. Now go to **8.4 X-ray Exposure** to start the exposure.

### 8.3.4 Waters' view Mode

#### Patient Positioning

1. Turn the Nasal Positioner to the **PA / Waters' view / Carpus** mode Positioning Marker as shown below.




2. Fold the Nasal Positioner up. The Nasal Positioner is not used in Waters' view mode.
3. Guide the patient to the CEPH unit.
4. Ask the patient to stand upright facing the sensor. Make sure that the patient's shoulders are level and that his/her neck is relaxed.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.



#### **WARNING**

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, correctly align the Ear Rods to the patient's ears, so his/her head does not move.
7. Direct the patient to swallow first before closing his/her mouth and guide the patient to bend the head backward  $30^{\circ}$  -  $40^{\circ}$ . Direct the patient to remain in the current position until image acquisition is completed.



8. Click the **READY** button. The x-ray exposure has not started yet.
9. Now go to **8.4 X-ray Exposure** to start the exposure.

### 8.3.5 Carpus Mode


For Carpus Mode, install the Carpus Plate first before positioning the patient.

#### Installing the Carpus Plate

1. Turn the Nasal Positioner to the **PA / Waters' view / Carpus** mode Positioning Marker as shown below.




2. Fold the Nasal Positioner up. The Nasal Positioner is not used in Carpus mode.
3. Fit the two ends of the Carpus Plate into the two holes of the CEPH unit as below.



4. Confirm that the Carpus Plate is safely mounted.

### Patient Positioning

1. Let the patient put his/her right hand splayed on the Carpus Plate as shown below. Make sure that the patient does not bend his/her fingers.



2. Ask the patient to close his/her eyes and stand still until the image acquisition is completed.
3. Click the **READY** button. The x-ray exposure has not started yet.
4. Now go to **8.4 X-ray Exposure** to start the exposure.

## 8.4 X-ray Exposure



- If an emergency occurs during image acquisition, release the **Exposure Switch** to cease X-ray emission.
- The operator shall always observe the X-ray safety regulations applicable to his/her area during the operation of this equipment.




- The operator must always keep vocal/visual contact with the patient during the image acquisition process.
- Do not operate the PC during exposure. Doing so may cause the system to malfunction.

**IMPORTANT**

- Let the patient close the eyes during the operation.
- To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also, don't let the patient move until the Temple Supports are open.

1. Get out of the X-ray room and close the door.
2. Press and hold down the **Exposure Switch** until image acquisition is completed.



• Yellow : X-ray On

**NOTICE**

The image appears on the screen.

**NOTICE**

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns yellow and “X-ray” changes to “X-RAY ON.”



X-RAY ON

3. Release the **Exposure Switch** when “Image capturing is completed” message appears on the screen.

## 8.5 Finishing the Scan

1. Leave enough space between the Ear Rods.
2. Fold the Nasal Positioner up in case it's unfolded.
3. Guide the patient out of the equipment.

## 8.6 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i / EasyDent**.

### **NOTICE**

Refer to the **EzDent-i / EasyDent User Manual** for more information.

1. The images are transferred to **EzDent-i / EasyDent** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.

# 9. Acquiring CBCT Images

## 9.1 CBCT Imaging Program Overview

- **Result Images**



It provides conventional 3D CT sliced images.



- **Image Acquisition Method**



It acquires images with the X-ray beam scanning specific oral & maxillofacial regions and reconstructs them to 3D sliced images.


- **Examination Programs**

It is classified as below based on the FOV.

Available FOVs (cm)	ROI	Description
16x9	 <p data-bbox="591 1041 797 1064">&lt;Occlusion / Center&gt;</p>	<ul style="list-style-type: none"> <li>- Covers full arch region, sinus and left/right TMJ.</li> <li>- Suitable for most oral surgery cases as well as multiple implant surgery.</li> </ul>
12x9	 <p data-bbox="591 1277 797 1300">&lt;Occlusion / Center&gt;</p> <p data-bbox="488 1499 625 1522">&lt;TMJ / Right&gt;</p> <p data-bbox="749 1499 886 1522">&lt;TMJ / Left&gt;</p> <p data-bbox="611 1711 776 1734">&lt;Airway / Center&gt;</p>	<ul style="list-style-type: none"> <li>- Covers both maxillary and mandibular structures including the 3<sup>rd</sup> molar region.</li> <li>- TMJ Right/Left and Airway mode are available.</li> </ul>

<p>8x9</p>	 <p>&lt;Occlusion / Right&gt;      &lt;Occlusion / Left&gt;</p> <p>&lt;Occlusion / Center&gt;</p> <p>&lt;TMJ / Right&gt;      &lt;TMJ / Left&gt;</p>	<ul style="list-style-type: none"> <li>- Selectable FOV for region (left/center/right).</li> <li>- Covers both maxillary and mandibular areas and left/right TMJ.</li> </ul>
<p>8x5</p>	 <p>&lt;Mandible / Center&gt;</p> <p>&lt;Mandible / Left&gt;      &lt;Mandible / Right&gt;</p> <p>&lt;Maxilla / Center&gt;</p>	<ul style="list-style-type: none"> <li>- Selectable FOV for region (left/center/right).</li> <li>- Covers both maxillary and mandibular areas.</li> </ul>

	 <p>&lt; Maxilla / Left&gt;      &lt; Maxilla / Right&gt;</p>	
<p>5x5</p>	 <p>&lt;Maxilla / Right Molar&gt; &lt;Maxilla / Left Molar&gt;          &lt;Maxilla / Right&gt; &lt;Maxilla / Left&gt;          &lt;Maxilla / Incisor&gt; &lt;Mandible / Incisor&gt;          &lt;Mandible / Right Molar&gt; &lt;Mandible / Left Molar&gt;          &lt;Mandible / Right&gt; &lt;Mandible / Left&gt;</p>	<ul style="list-style-type: none"> <li>- Covers 3 ~ 4 areas through capturing ROI.</li> <li>- Can acquire 3 ~ 4 teeth at once.</li> <li>- Endo mode (Voxel Size: 0.08) applied.</li> </ul>


<p>Double Scan (Optional)</p>	 <p>1. &lt;Mandible / Center&gt;</p> <p>2. &lt;Maxilla / Center&gt;</p>	<ul style="list-style-type: none"> <li>- Full arch, sinus, left and right TMJ area can be checked</li> <li>- Suitable for most intra-oral surgery including multiple implant placement</li> </ul>
	<div style="border: 1px dashed black; padding: 5px;"> <p><b>NOTICE</b></p> <p>After taking 2 consecutive exposures in the order of "1 → 2", it is automatically synthesized and displayed as one full image.</p> </div>	

## 9.2 Configuring Exposure Parameters

To acquire CBCT Images, **6. Getting Started** must be completed first.

### NOTICE

You can set the imaging parameters on either the Touchscreen or the Console Software running on the PC. They are synchronized and display the same environmental settings.





- ◆ Available options for each FOV are as below.


Available FOV (cm)	Vertical option	Horizontal option	Tooth option
16x9	Occlusion	Center	N/A
12x9	Occlusion	Center	N/A
	TMJ	Right	
		Left	
Airway	Center		
8x9	Occlusion	Right	N/A
		Center	
		Left	
	TMJ	Right	
Left			
8x5	Maxilla	Right	N/A
		Center	
		Left	
	Mandible	Right	N/A
		Center	
Left			
5x5	Maxilla / Mandible	N/A	Right Molar
			Right
			Incisor
			Left
			Left Molar
16x15 (Double Scan; 16x9 stitching) (Optional)	Mandible (First Scan) Maxilla (Second Scan)	Center	N/A

3. Select an Image Option. (On Touch Screen, click **Settings** button before selecting options.)



4. Select a Voxel Size.

**NOTICE** MAR (Metal Artifact Reduction) function is applied automatically if there are metal objects in the image. MAR may increase image reconstruction time.



5. The Gender / Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.



**NOTICE**

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

6. Select X-ray intensity.



<Console Software>



<Touch Screen>

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft:

Soft ≤ Normal ≤ Hard


**NOTICE**

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Inten
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

7. The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ±1 kVp and ± 0.1 mA respectively.




<Console Software>




<Touch Screen>

8. Select the Process Option. The HD is the default image reconstruction process option of the equipment. The Fast is noisier than the HD, but the CT image is reconstructed quickly. (With MAR, display metal better.)



<Console Software>



<Touch Screen>

9. Click the **CONFIRM** button when the exposure parameter setting is completed.



<Console Software>



<Touch Screen>

**NOTICE**

When you click **CONFIRM** button,

- The Rotating Unit will move to its initial scanning position.
- The Vertical Beam will be activated to make patient positioning easier.
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below the Patient Information window.

DAP  
127.334307 mGy x cm<sup>2</sup>

Scan-time  
13.5 Sec

Exposure-time  
13.5 Sec

10. Guide the patient to the equipment.

### 9.3 Obtaining Double Scan Image (optional)


FOV 16 × 9 images can be acquired by upper/lower double scan and stitched to acquire full FOV 16x15 CBCT images.

#### 9.3.1 First Scan Step of Double Scan

To get a double scan image, first **6. Getting Started** should be completed.

### NOTICE

The PC's Console Software has the same function as the touch screen of the device. The Console Software and the touch screen are linked to each other in real-time, and the environment settings are the same.



1. Click the Double Scan button on the CBCT main screen.




<Console Software>




<Touchscreen>

2. Select the Image Option. (On the touchscreen, first click the Settings button.)



< Console Software >



< Touchscreen >

3. The patient's gender/age group is automatically selected based on patient information. You can choose manually if necessary.



<Console Software>



<Touchscreen>

**NOTICE**

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

4. Select the X-ray intensity.



<Console Software>



<Touchscreen>

X-ray intensity can be classified as Hard, Soft, Normal according to the patient's head circumference.

Soft ≤ Normal ≤ Hard

**NOTICE**

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

5. The tube voltage and tube current values are automatically set according to the patient's performance and X-ray intensity. If you click the arrow button, you can make fine adjustments in ± 1 kVp, ± 0.1 mA increments.



<Console Software>



<Touchscreen>

- After setting the Exposure conditions, click the **CONFIRM** button.



<Console Software>



<Touchscreen>

**NOTICE**


When you click the **CONFIRM** button

- The rotating unit moves to the initial position.
- The vertical beam is activated.
- Below the patient information window, DAP (Dose Area Product), Scan Time, Exposure Time values are displayed.

DAP	1435.017724 mGy x cm <sup>2</sup>
Scan-time	9.0 Sec
Exposure-time	9.0 Sec

- Guide the patient to the equipment and align the patient's posture. For more information on patient posture alignment, see section **9.3 Patient Positioning**
- Press the **READY** button. At this time, X-rays are not irradiated.
- Go to **9.5 X-ray Exposure** and proceed with X-ray exposure.
- Guide the patient out of the shielded room with equipment.

### 9.3.2 Second Scan Steps of Double Scan



1. Select the **Maxilla** as the **Vertical Option**. It is automatically selected.



2. The tube voltage and tube current values are automatically set according to the patient's performance and X-ray intensity. If you click the arrow button, you can make fine adjustments in  $\pm 1$  kVp,  $\pm 0.1$  mA increments



<Console Software>

<Touchscreen>

3. Click the **CONFIRM** button



<Console Software>



<Touchscreen>

**NOTICE**

When you click the **CONFIRM** button

- The rotating unit moves to the initial position.
- The vertical beam is activated.

Below the patient information window, DAP (Dose Area Product), Scan Time, Exposure Time values are displayed.

DAP  
1435.017724 mGy x cm<sup>2</sup>

Scan-time  
9.0 Sec

Exposure-time  
9.0 Sec

4. Guide the patient to the equipment and align the patient's posture. For more information on patient posture alignment, see section **9.3 Patient Positioning.**
5. Press the **READY** button. At this time, X-rays are not irradiated.
6. Go to **9.5 X-ray Exposure** and proceed with X-ray exposure.
7. Guide the patient out of the shielded room with equipment.
8. Move to the **9.6 Finishing the Scan** to finish shooting and check the image.

## 9.4 Patient Positioning




- Have patience (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation.
- Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

### IMPORTANT

- Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition.
- Metal implants or bridges may reduce the quality of the images.
- Be sure to adjust the laser beam correctly. Otherwise, the quality of images can be lower due to ghost images or expansion/reduction of the images.

### Getting prepared


1. Let the patient remove all the metal objects (glasses, earrings, hairpins, braces, false teeth, etc.). Metal objects may induce ghost images and lower image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button or switch option to adjust the equipment to match the height of the patient.



< Control Panel – LCD type >

### Normal Patient Positioning

1. Insert the Normal Bite into the Normal Chinrest and cover it with a Sanitary Vinyl Cover.



- The Sanitary Vinyl Cover is for single use only. It should be replaced for each patient. Be sure to use the approved vinyl cover.
- 
- Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.


2. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.



3. Guide the patient to the equipment.
4. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.
5. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.
6. Let the patient bite the Bite along its grooves with his/her front teeth.




7. Let the patient maintain the posture as follows:
- Close the mouth.
  - Place the tongue to the roof of the mouth.
  - Close the eyes.



### Double Scan (Mandible) Patient Positioning


1. Remove the **Special Chinrest** and insert the **Normal Chinrest** into the equipment.
2. Insert the **Special Bite B** into the **Normal Chinrest**.



3. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.



4. Guide the patient to the equipment.
5. Use the **Vertical Frame Up/Down** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.



6. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.
7. Let the patient maintain the posture as follows:
  - Close the eyes.


#### **NOTICE**

- Ask the patient to maintain his/her position until the operation is completed.

### **Double Scan (Maxilla) Patient Positioning**

1. Remove the **Normal Chinrest** and **Special Bite B**.
2. Insert the **Double Scan Support** into the equipment then cover with a sanitary vinyl cover.

Double Scan Support



**WARNING**




The sanitary vinyl cover is for single use only. It should be replaced after each patient. Be sure to use the approved vinyl cover.

3. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.



4. Guide the patient to the equipment.
5. Use the **Vertical Frame Up/Down** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.



6. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.


7. Let the patient maintain the posture as follows:
  - Close the eyes.

**NOTICE**

- Ask the patient to maintain his/her position until the operation is completed.

**TMJ Patient Positioning**

1. Remove the **Normal Chinrest** and insert the **Special Chinrest** into the equipment.
2. Insert the **Special Bite A** into the **Special Chinrest**



3. Use the **Temple Supports OPEN/CLOSE** button on the control panel to widen the Temple Supports.



4. Guide the patient to the equipment.
5. Use the **Vertical Frame Up/Down** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.
6. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
  - Hold the handles tightly.
  - Press the chest against the equipment.
  - Keep both feet close to the inside of the base.
  - Keep both shoulders parallel.
  - Straighten the Cervical Spine and stand still.
7. Let the patient maintain the posture as follows:
  - Close the eyes.

**NOTICE**

- Ask the patient to maintain his/her position until the operation is completed.

### Laser Beam Aligning




Be careful not to shine the Laser Beam directly into the person's eyes. Doing so may result in vision loss.



If the Laser Beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align Laser Beam properly.

Vertical Beam  
/ Mid-sagittal plane

Horizontal Beam  
/ Frankfurt plane



e.g. FOV 16 X 9

### NOTICE

This is a sample illustration for reference only. Actual FOV may vary from the image as shown above.

### IMPORTANT

- The horizontal beam is used to point to the center of the height within the FOV.
- Unlike PANO mode, the horizontal beam does not move only in the CT mode. this is not a malfunction/error.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It's to prevent the horizontal expansion of the image)

### Finishing Patient Positioning

1. After checking the positions of the patient and the Laser Beam, click the **Temple Supports OPEN/CLOSE** button on the control panel to prevent the patient's head from moving.



Make sure that the Temple Supports are in the CLOSE position before clicking the **READY** button.

2. Click the **READY** button. X-ray exposure has not started yet.
3. Now go to **9.5 X-ray Exposure** to start the exposure.

## 9.5 X-ray Exposure



- If an emergency occurs during image acquisition, release the **Exposure Switch** to cease X-ray emission.
- The operator shall always observe the X-ray safety regulations applicable to his/her area during the operation of this equipment.



- The operator must always keep vocal/visual contact with the patient during the image acquisition process.
- Do not operate the PC during exposure. Doing so may cause the system to malfunction.

**IMPORTANT**


- Let the patient close the eyes during the operation.
- To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also, don't let the patient move until the Temple Supports are open.

1. Get out of the X-ray room and close the door.

**IMPORTANT**

The operator must always keep vocal/visual contact with the patient during image acquisition.

2. Press and hold down the **Exposure Switch** until image acquisition is completed.



• Yellow : X-ray On

**NOTICE**

The image appears on the screen.

**NOTICE**

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns yellow and "X-ray" changes to "X-RAY ON."



X-RAY ON

3. Release the **Exposure Switch** when "Image capturing is completed" message appears on the screen.

## 9.6 Finishing the Scan

1. Open the Temple Supports and guide the patient out of the equipment.
2. For Normal Bite, remove the Sanitary Vinyl Cover from the Bite.
3. Press the **READY** button to bring the Rotating Unit back to its initial position.

## 9.7 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i / EasyDent**.

### NOTICE

Refer to the **EzDent-i / EasyDent User Manual** for more information.

1. The images are transferred to **EzDent-i / EasyDent** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.
4. Then, **Ez3D-i / Ez3D Plus** will run automatically for 3D viewing.

Left blank intentionally

## 10. Acquiring 3D MODEL Scan Images

### 10.1 3D MODEL Scan Imaging Program Overview

- **Result Images**


It provides 3D modeling surface data of the Plaster Cast. (STL file)

- **Image Acquisition Method**

It acquires images with the X-ray beam scanning the Plaster Cast and reconstructs them to 3D sliced images and converts the sliced images into 3D modeling surface data.

- **Examination Programs**

It is classified as below based on the MODEL type.

Applied FOV (cm)	Vertical Option	ROI	Description
8x9	Upper (Maxilla)		Captures a whole maxillary Plaster Cast.
	Lower (Mandible)		Captures a whole mandibular Plaster Cast.

#### **NOTICE**


The 3D MODEL Scan modality is not available for EasyDent / Ez3D Plus users.

## 10.2 Configuring Exposure Parameters


To acquire 3D MODEL Scan Images, **6. Getting Started** must be completed first.

### NOTICE

You can set the imaging parameters on either the Touch Screen or the Console Software running on the PC. They are synchronized and display the same environmental settings.



1. Click the MODEL button on the Main Screen.



<Console Software>



<Touch Screen>


**NOTICE**

The **CEPH** button exists only when the CEPH imaging program is included in the equipment.

2. Select Model Examination type.



<Console Software>



<Touch Screen>

3. The Gender / Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.



<Console Software>



<Touch Screen>

**NOTICE**

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

4. Select X-ray intensity.



<Console Software>



<Touch Screen>

5. The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by  $\pm 1$  kVp and  $\pm 0.1$  mA respectively.



<Console Software>



<Touch Screen>

6. Click the **CONFIRM** button when the exposure parameter setting is completed.



<Console Software>




<Touch Screen>

**NOTICE**

When you click **CONFIRM** button,

- The Rotating Unit will move to its initial scanning position.
- The Vertical Beam will be activated to make patient positioning easier.
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below the Patient Information window.




7. Bring the Plaster Cast to the equipment.

## 10.3 MODEL Positioning

### MODEL Scan Jig Installation

1. Remove the Temple Supports and the Chinrest
2. Insert the MODEL Scan Jig.



### Laser Beam Aligning

1. Put the Plaster Cast on the MODEL Scan Jig. (Whether the Plaster Cast is for Maxilla or Mandibular, place it flat side down.)



2. Align the Mid-sagittal plane Laser Beam to the center of the Plaster Cast. (To prevent the horizontal expansion of the image)




3. Click the **READY** button. X-ray exposure has not started yet.
4. Now go to **10.4 X-ray Exposure** to start the exposure.

## 10.4 X-ray Exposure



Do not operate the PC during exposure. Doing so may cause the system to malfunction.

1. Get out of the X-ray room and close the door.
2. Press and hold down the **Exposure Switch** until image acquisition is completed.



• Yellow : X-ray On

### NOTICE

The image appears on the screen.

### NOTICE

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns yellow and "X-ray" changes to "X-RAY ON."



X-RAY ON

3. Release the **Exposure Switch** when "Image capturing is completed" message appears on the screen.
4. Remove the Plaster Cast out of the equipment.

## 10.5 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM or STL (Stereo Lithography) format.

**NOTICE**

Refer to the **EzDent-i User Manual** for more information.

1. The images are transferred to **EzDent-i** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.
4. You can check the captured image with a 3<sup>rd</sup> party STL viewer.

Left blank intentionally

# 11. Troubleshooting

## 11.1 Troubleshooting

If a problem occurs while operating the equipment, perform the corresponding troubleshooting measures outlined in the table below. If the problem persists, please contact our customer support staff.

### If the equipment is not working

Cause	Actions to be taken
Failure of power supply	Check the equipment's power supply.
Initialization status	Wait until the equipment has been initialized and then try again.
Failure of the Control PC's connection	Check the connection status of Communication Port (Optic) which connects the PC to the equipment.

### If the Exposure Switch is not functioning

Cause	Actions to be taken
Failure of readiness	Check whether the Console Software is ready for imaging.

### If imaging cannot be performed

Cause	Actions to be taken
Failure of initialization	Wait until the equipment is initialized and then try again. If this problem persists, restart the equipment.

### If the Laser Beam has shut off and patient positioning cannot be performed

Cause	Actions to be taken
Expiration of the time allotted for patient positioning	Press the Laser Beam button to turn on the Laser Beam.


## 11.2 Error Codes

In instances of abnormal operation, error messages will be displayed with error codes on the Console Software and Control Panel. If a problem persists, please request assistance from the customer support information services.

Error messages will be displayed in the format written below.

**[Code: E00X.0XX.0XX]**

The code consists of three parts: Main code, Subcode, Tail code.



### NOTICE

- **The main code** indicates the source of error codes. The source is categorized as hardware, software, an acquisition module, etc.
- **Subcode** describes the specific area where the error has occurred according to the main code.
- **The tail code** explains the specific symptoms and causes of the errors mentioned in the subcode.

### 11.2.1 Main code - Hardware (001)

#### 11.2.1.1 Subcode – Generator related error (001)

Tail code	Description	Solution
001	Appears when the tube is not ready for use	1) Check the CAN communication response by sending the following command to the inverter board, sensor, and the collimator to check the CAN communication operation of the Main MCU. <ul style="list-style-type: none"> <li>- Inverter board: [SPM_IVER]</li> <li>- Sensor: [SPM_FISS_0001]</li> <li>- Collimator: [SPM_CVER]</li> </ul> 2) If the inverter board, sensor, and the collimator do not respond, replace the Main MCU. 3) If only the inverter board does not respond, check the connection status of H001909A (Pin 1, 3) and H001922A (Pin 7, 8) cable to check CAN communication between the Main MCU and inverter board.

Tail code	Description	Solution
		<p>※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>3-1) If the cable is normal, check that 24V is normally applied to pins 1, 2 of CN13 to check the input voltage of the inverter board.</p> <p>3-2) When 24V is normally applied to the inverter board, replace the generator.</p> <p>3-3) If 24V is not confirmed, check the power cable H001909A (Pin 21, 22) and H001922A (Pin 1, 2) from the Main MCU to the inverter board.</p> <p>3-4) If the cable is normal, see if DC 24V is normally applied to the pins 4, 5 of CN1600 to check the input power of the Main MCU.</p> <p>3-5) If 24V is applied, replace the Main MCU.</p> <p>3-6) If 24V is not applied, replace the power board.</p>
002	Appears when cable between tube tank and Inverter board are disconnected	<p>1) Check the connection of high voltage cable between the inverter board and the tube tank.</p> <p>2) If the same problem occurs after cable replacement, replace the generator.</p> <p>※ Generator includes a tube tank and the inverter board.</p> <p>※ Generator inverter high voltage cable: H001907A</p>
003	Appears when a current of the inverter board exceeds maximum allowable level during X-ray irradiation	<p>1) Send out the following command to the motor, inverter board, sensor, and the collimator to check the operation status and IC board component status. (Advance remote check required)</p> <ul style="list-style-type: none"> <li>- Motor: [SPM_TEST]</li> <li>- Inverter board: [SPM_IIVER]</li> <li>- Sensor: [SPM_FISS_0001]</li> <li>- Collimator: [SPM_CVER]</li> </ul> <p>- If there is a problem during operation or communication, replace the corresponding component (motor, board) and generator.</p> <p>- If there is no problem, replace the generator.</p> <p>2) Check the high voltage input power (power board</p>

Tail code	Description	Solution
		<p>output voltage) of the inverter board CN11 pin 1, 2 is DC 380V±10%.</p> <p>3) When X-ray is irradiated, check whether the high voltage input power of the inverter board CN11 pin 1, 2 is 340V or more. (Normal if 340V or more)                      - If any problem occurs in stage 2 or 3, replace the power board.</p> <p>4) For equipment manufactured prior to 2014/04/19, apply the brass bar solution.</p>
004	Appears when there is ±10kV or more voltage difference in tube voltage compared to a reference value	Follow the instructions described in the tail code 003.
005	Appears when there is ±0.5mA or more current difference in tube current compared to a reference value	Follow the instructions described in the tail code 003.
006	Appears when there is ±20kV or more voltage difference in tube voltage feedback compared to the average value	Follow the instructions described in the tail code 003.
007	Appears when there is ±1mA or more current difference in tube current feedback compared to the average value	Follow the instructions described in the tail code 003.
008	Appears when the temperature of the mono tank is	1) Check that the temperature of mono-tank is 55°C or higher through the log message of the capture

Tail code	Description	Solution
	above the setting temperature	<p>program.</p> <p>※ The temperature is automatically checked every 10 seconds in the terminal program.</p> <p>2) If the temperature is higher than 55°C, the problem has occurred from continuous shooting using the capture program. It is necessary to notify and educate the staff members that enough cooling time is required during shooting.</p> <p>※ The cooling time Vatech suggests is 1:60. (60 seconds pause at the time of 1-second irradiation)</p> <p>3) If the same problem occurs after the measurements stated above, replace the generator.</p>
009	Appears when the inverter output current is higher than 1A during X-ray irradiation (In EP, IP condition)	Follow the instructions described in the tail code 003.
010	Appears when inverter board falsely recognizes the exposure switch signal as OFF after the irradiation On command	<p>1) Check if 24V is typically applied to the inverter board.</p> <p>- If there is any problem with the power, check if 24V is typically outputted from H001922A (Main MCU), pins 1, 2 of CN13 (inverter board), and pins 4, 5 of CN3 (power board).</p> <p>2) Check the multiple cables. Check the connection status of H001913A of the exposure switch from the Main MCU and the inverter board.</p> <p>※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>3) If the problem is not resolved, replace the generator.</p>
011	Appears when X-	1) Make sure that the staff member correctly presses

Tail code	Description	Solution
	ray OFF command is not sent to inverter board in 0.5 seconds after turning off the exposure switch	the exposure switch during irradiation. 2) Replace the exposure switch and check if the same error occurs. 3) Check the connection status of the H001909A (pin 1, 3) and H001922A (pin 7, 8) from the Main MCU to the inverter board. ※ Connection status refers to reconnection, disconnection, pin status, etc. 3-1) If the problem is not resolved, replace the multiple cables. 4) Check if 24V is applied to the Main MCU and the inverter board. 4-1) Check if the input voltage of pins 5, 6 of CN600 of the Main MCU is 24V. - If not 24V, replace the power board. 4-2) Check if the input voltage of pin 1, 2 of CN13 of the inverter board is 24V. - If not 24V, replace the Main MCU. - If 24V is applied regularly, replace the generator.
012	Appears when kV feedback is over -20kV compared to the setting value during X-ray irradiation	Follow the instructions described in the tail code 003.
013	Appears when kV feedback is over +20kV compared to the setting value during X-ray irradiation.	Follow the instructions described in the tail code 003.
014	Appears when the mA feedback value is less than 50% compared to setting conditions during X-ray irradiation.	Follow the instructions described in the tail code 003.
015	Appears when the	Follow the instructions described in the tail code 003.

Tail code	Description	Solution
	mA feedback value is higher than 150% compared to setting conditions during X-ray irradiation.	

### 11.2.1.2 Sub code - Motor related error (002)

Tail code	Description	Solution
020	Appears during p-axis motor origin movement	<p>1) Check the operation status due to mechanical load.            ※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>2) Check if the photosensor LED of the p-axis motor and the XORG LED of the Main MCU lights up. Also, check if 24V is outputted from pins 4, 5 of CN1701.            - If there is a problem, replace the photosensor.</p> <p>3) If p-axis motor usually operates, increase the torque value of the p-axis motor.            ※ Back up the torque value by [SPM_XMRT] and apply the torque value as 2000 higher than the current value.</p> <p>4) Check the cable connection of the p-axis motor.            ※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>5) Check the p-axis motor for defects and replace the motor if a problem occurs.            5-1) Make sure the resistance between pins 1, 3 and pins 2, 4 of CN1701 is 1.1~1.3 Ω.            5-2) Check if there is a short between pins 1, 4 and pins 2, 3 of CN1701.</p> <p>6) If the problem is not resolved, replace the Main MCU.</p>
021	Appears during rotator-axis motor origin movement	<p>1) Check the operation status due to mechanical load.            ※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p>

Tail code	Description	Solution
		<p>2) Check if the photosensor LED of the rotator-axis motor and the RORG LED of the Main MCU lights up. Also, check if 24V is outputted from the pins 1, 2 of CN2708 (rotator) / CN1600 (Main MCU). - If there is a problem, replace the photosensor.</p> <p>3) If rotator-axis motor usually operates, increase the torque value of the rotator-axis motor. ※ Back up the torque value by [SPM_RMRT] and apply the torque value as 2000 higher than the current value.</p> <p>4) Check the cable connection of the rotator-axis motor. ※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>5) Check the rotator-axis motor for defects and replace the motor if a problem occurs.</p> <p>5-1) Make sure the resistance between pins 1, 3 and pins 2, 4 is 1.1~1.3 Ω. 5-2) Check if there is a short between pins 1, 4 and pins 2, 3.</p> <p>6) If the problem is not resolved, replace the Main MCU.</p>
027	Appears during Cephalo sensor motor origin movement	<p>1) Check the operation status due to mechanical load. ※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc. -If there is a problem, replace the photosensor.</p> <p>2) Check if the LED of the cephalo sensor motor lights up. Also, check if 24V is outputted from pins 1, 3 of CN302. - If there is a problem, replace the photosensor.</p> <p>3) Increase the torque value of the cephalo sensor motor. ※ Back up the torque value by [SPM_HMRT] and apply the torque value as 2000 higher than the current value.</p> <p>4) Check the cable connection of the cephalo sensor axis motor. ※ Connection status refers to reconnection, disconnection, pin status, etc.</p>

Tail code	Description	Solution
		<p>5) Check the cephalo sensor axis motor for defects and replace the motor if a problem occurs.</p> <p>5-1) Make sure that the resistance between pins 1, 3 and pins 2, 4 is 4.5~6Ω.</p> <p>5-2) Check if there is a short between pins 1, 4 and pins 2, 3.</p> <p>6) If the problem is not resolved, replace the main MCU.</p>
030	Appears during quadruple axis collimator Left origin movement	<p>1) Check if the collimator MCU FW is the latest version.</p> <p>※ When upgrading, the firmware upgrade cable is required.</p> <p>2) If the problem is not resolved, check the operation status due to mechanical load.</p> <p>※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>- If any foreign substance is found, remove the substance and lubricate the component.</p> <p>3) Apply the torque value as 100~200 higher than the current value.</p> <p>※ Check and back up the current torque value by [SPM_TQB?]. Adjust the torque value by [SPM_TOQB_xxxx].</p> <p>4) If the problem is not resolved, replace the collimator assembly.</p> <p>5) If the problem is not resolved after replacing the collimator assembly, replace the Main MCU.</p>
031	Appears during quadruple axis collimator Right origin movement	<p>1) Check if the collimator MCU FW is the latest version.</p> <p>※ When upgrading, the firmware upgrade cable is required.</p> <p>2) If the problem is not resolved, check the operation status due to mechanical load.</p> <p>※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical</p>

Tail code	Description	Solution
		<p>load, timing belt, etc.</p> <p>- If any foreign substance is found, remove the substance and lubricate the component.</p> <p>3) Apply the torque value as 100~200 higher than the current value.</p> <p>※ Check and back up the current torque value by [SPM_TQA?]. Adjust the torque value by [SPM_TOQA_xxxx].</p> <p>4) If the problem is not resolved, replace the collimator assembly.</p> <p>5) If the problem is not resolved after replacing the collimator assembly, replace the Main MCU.</p>
032	Appears during quadruple axis collimator Up origin movement	<p>1) Check if the collimator MCU FW is the latest version.</p> <p>※ When upgrading, the firmware upgrade cable is required.</p> <p>2) If the problem is not resolved, check the operation status due to mechanical load.</p> <p>※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>- If any foreign substance is found, remove the substance and lubricate the component.</p> <p>3) Apply the torque value as 100~200 higher than the current value.</p> <p>※ Check and back up the current torque value by [SPM_TQD?]. Adjust the torque value by [SPM_TOQD_xxxx].</p> <p>4) If the problem is not resolved, replace the collimator assembly.</p> <p>5) If the problem is not resolved after replacing the collimator assembly, replace the Main MCU.</p>
033	Appears	1) Check if the collimator MCU FW is the latest version.

Tail code	Description	Solution
	during quadruple axis collimator down origin movement	<p>※ When upgrading, the firmware upgrade cable is required.</p> <p>2) If the problem is not resolved, check the operation status due to mechanical load.</p> <p>※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>- If any foreign substance is found, remove the substance and lubricate the component.</p> <p>3) Apply the torque value as 100~200 higher than the current value.</p> <p>※ Check and back up the current torque value by [SPM_TQC?]. Adjust the torque value by [SPM_TOQC_xxxx].</p> <p>4) If the problem is not resolved, replace the collimator assembly.</p> <p>5) If the problem is not resolved after replacing the collimator assembly, replace the Main MCU.</p>
037	Appears during generator tilting	<p>1) Make sure the product is a Ceph model. If it is not a Ceph model, check whether [SPM_CISC], [SPM_TITY] is activated.</p> <p>2) If the output is standard, check the operation status due to mechanical load.</p> <p>※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>- After sending the [SPM_TICE] command, check if the generator is tilting in Ceph mode and the limit switch is detected. After confirmation, send [SPM_TIFR] to stop the motor drive.</p> <p>- After sending [SPM_TIPA] command, check whether the generator is tilting in the PANO/CT mode and the limit switch is detected. After confirmation, send [SPM_TIFR] to stop the motor drive.</p> <p>3) If the limit switch is not recognized usually, replace the</p>

Tail code	Description	Solution
		<p>Main MCU.</p> <p>4) If the tilting motor does not operate normally, check the tilting motor for defects. 4-1) If the tiling motor is standard, please replace the Main MCU.</p>
038	Appears during temple support the motor operation	<p>1) If there is any noise during operation, replace the ass'y.</p> <p>2) Check the limit switch of the temple support. ※ Check whether the TSO1 and TSO2 LED light up.</p> <p>3) Check the ORG cable of the temple support.</p> <p>4) Check the alignment of the motor and the bearing block. (critical*)</p> <p>5) Adjust the torque value and see if the problem is resolved. ※ Check the current torque value by [SPM_TQK?]. Adjust the torque value by [SPM_TOQK_xxxx].</p> <p>6) If the same problem occurs, replace the Main MCU.</p> <p>7) If the problem is not resolved, replace the temple support ass'y.</p>
039	Appears during X-axis motor origin movement	<p>1) Check the operation status due to mechanical load. ※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>2) Check if the LED of the x-axis motor and the XORG LED of the Main MCU lights up. - If there is a problem, replace the photosensor.</p> <p>3) If the x-axis motor usually operates, increase the torque value of the x-axis motor. ※ Back up the [SPM_XMRT] value and apply 1 increment of torque.</p>

Tail code	Description	Solution
		<p>4) Check the cable connection of the x-axis motor.            ※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>5) Check the x-axis motor for defects and replace the motor if a problem occurs.            5-1) Make sure the resistance between pins 1, 3 and pins 2, 4 is 1.1~1.3Ω.            5-2) Check if there is a short between pins 1, 4 and pins 2, 3.</p> <p>6) If the problem is not resolved, replace the Main MCU.</p>
040	Appears during Y-axis motor origin movement	<p>1) Check the operation status due to mechanical load.            ※ Operation status includes motor operation, a foreign substance in the axis, the existence of the mechanical load, timing belt, etc.</p> <p>2) Check if the LED of the y-axis motor and the YORG LED of the Main MCU lights up.            - If there is a problem, replace the photosensor.</p> <p>3) If the y-axis motor usually operates, increase the torque value of the Y-axis motor.            ※ Back up the [SPM_YMRT] value and apply 1 increment of torque.</p> <p>4) Check the cable connection of the Y-axis motor.            ※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>5) Check the y-axis motor for defects and replace the motor if a problem occurs.            5-1) Make sure that the resistance between pins 1, 3 and pins 2, 4 is 1.1~1.3Ω.            5-2) Check if there is a short between pins 1, 4 and pins 2, 3.</p> <p>6) If the problem is not resolved, replace the Main MCU.</p>

## 11.2.1.3 Subcode – Exposure switch related error (003)

Tail code	Description	Solution
060	It appears if the exposure switch is pressed when turning on the equipment.	<p>1) To check whether the exposure switch usually operates or not, conduct the following measurement.</p> <p>1-1) Remove the exposure switch and see if the same problem occurs. 1-2) Replace the current exposure switch with an extra switch provided - Replace the exposure switch if the problem occurs after conducting the measurements mentioned above.</p> <p>2) Check the status of the extension cable. ※ Connection status refers to short-circuited, disconnection, etc.</p> <p>3) If there is no problem in stages 1 and 2, replace the Main MCU.</p>

## 11.2.1.4 Subcode – Other error (004)

Tail code	Description	Solution
102	Appears when there is no response during CAN communication.	<p>1) Replace the Main MCU and the collimator manufactured prior to 2019.01.17, S/N: 069-002730 when H102 error occurs intermittently. If an error occurs all the time or after replacing the Main MCU and the collimator, follow the instructions written below.</p> <p>2) Check the CAN communication response by sending the following command to the inverter board and the sensor to check the CAN communication operation of the Main MCU. - Inverter board: [SPM_IVER] - Sensor: [SPM_FISS_0001]</p> <p>3) If both the inverter board and sensor do not respond, replace the Main MCU.</p> <p>4) If only the inverter board does not respond, check the connection status of CN13 (H001922A cable of the</p>

Tail code	Description	Solution
		<p>inverter) to check CAN communication between the Main MCU and the inverter board.</p> <p>※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>4-1) If the cable is normal, check that 24V is normally applied to pins 1, 2 of CN13 to check the input voltage of the inverter board.</p> <p>4-2) When the 24V power supply is normal, replace the inverter board.</p> <p>4-3) If 24V is not confirmed, check the H001904A cable from the Main MCU to the inverter board.</p> <p>4-4) If the cable is normal, check if 24V is applied to pins 1, 2 of CN13 to check the input power of the Main MCU.</p> <p>4-5) If 24V is applied, replace the Main MCU.</p> <p>4-6) If 24V is not applied, replace the power board.</p> <p>5) If only the sensor does not respond, check the connection status of H001909A, H001918A, and H001925A cable to check CAN provide communication between the main MCU and the sensor.</p> <p>※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>5-1) If the cable is normal, check that 8V is normally applied to the sensor.</p> <p>5-2) When the 8V power supply is normal, replace the sensor.</p> <p>5-3) If 8V power is not confirmed, check if 24V is normally applied to pins 6, 8 of CN2801 (DC04 board).</p> <p>5-4) If the 8V is applied, replace the DC04 board.</p> <p>5-5) If the 8V is not confirmed, check if 24V is normally applied to pins 21, 22 of CN2704 (SUB046 board).</p> <p>5-6) If the 24V is applied, replace the SUB046 board.</p> <p>5-7) If the 24V is not confirmed, check if 24V is normally applied to pins 5, 6 of CN1600 (Main MCU board).</p> <p>5-8) If 24V is applied, replace the Main MCU.</p> <p>5-9) If 24V is not confirmed, replace the power board.</p> <p>6) If only the collimator does not respond, check the connection status of H001909A and H001915A cable to check CAN communication between the Main MCU</p>

Tail code	Description	Solution
		<p>and the collimator.</p> <p>※ Connection status refers to reconnection, disconnection, pin status, etc.</p> <p>6-1) If the cable is normal, check if 24V is normally applied to pins 1, 4 of CN500.</p> <p>6-2) If 24V is applied, replace the collimator ass'y.</p> <p>6-3) If 24V is not confirmed, check if 24V is normally applied to pins 21, 22 of CN2704 (SUB046 board).</p> <p>6-4) If 24V is applied, replace the SUB046 board.</p> <p>6-5) If 24V is not confirmed, check if 24V is applied to pins 5 and 6 of CN1600 (Main MCU board).</p> <p>6-6) If 24V is applied, replace the Main MCU.</p> <p>6-7) If 24V is not confirmed, replace the power board.</p>

## 11.2.2 Main code – Software (002)

### 11.2.2.1 Sub code – Sequence related error (001)

Tail code	Description	Solution
001	Appears when the packing mode is enabled	Check whether the packing mode is enabled.
002	Appears when the door is open	<p>1) Check whether the door is open.</p> <p>2) Check if the activation of the DoorLock function in the capturing SW and the firmware of the main MCU.</p> <p>3) Check the [SPM_DROP_000x] value by sending [SPM_ISDR] through serial communication. (0: Deactivated, 1: Activated)</p> <p>4) Check the door sensor.</p> <p>5) If the problem is not resolved, contact the customer service team for further information.</p>
003	Appears when the exposure switch is pressed	<p>1) To check whether the exposure switch usually operates or not, conduct the following measurement.</p> <p>1-1) Remove the exposure switch and see if the same</p>

Tail code	Description	Solution
		<p>problem occurs.</p> <p>1-2) Replace the current exposure switch with an extra switch provided</p> <ul style="list-style-type: none"> <li>- Replace the exposure switch if the problem occurs after conducting the measurements mentioned above.</li> </ul> <p>2) Check the status of the extension cable.</p> <p>※ Connection status refers to short-circuited, disconnection, etc.</p> <p>3) If there is no problem in stages 1 and 2, replace the Main MCU.</p>

#### 11.2.2.2 Sub code – PC Resolution related error (010)

Tail code	Description	Solution
001	Appears when the resolution is less than 1280x1024	<p>1) Check whether the resolution setting is 1280x1024 or higher than 1280x1024.</p> <p>2) Check whether the magnification of the resolution is 100%.</p>
002	Appears when the resolution is less than 1200x960	<p>1) Check whether the resolution setting is 1200x960 or higher than 1200x1024.</p> <p>2) Check whether the magnification of the resolution is 100%.</p>

#### 11.2.2.3 Subcode – PC Network related error (024)

Tail code	Description	Solution
002	Appears when the port is invalid	<p>1) Check whether the power of the equipment is on when the capture SW is running.</p> <p>2) Check the communication port setting. (Default COM4)</p> <p>3) Check the communication cable connection status.</p> <p>※ Connection status refers to reconnection, bending, contamination, etc.</p>
003	Appears when the time is out	<p>1) Check whether the power of the equipment is on</p>

Tail code	Description	Solution
		<p>when the capture SW is running.</p> <p>2) Check the optic cable of each section from PC to the Main MCU by using the extra optic cable provided.</p> <p>3) If the cable is regular, replace the Main MCU.</p>

### 11.2.3 Main code - Acquisition Module (003)

#### 11.2.3.1 Sub code – Initialization Failure related error (010)

Tail code	Description	Solution
000	Appears when COM port cannot be opened	<p>1) Check whether the driver of the grabber and virtual comport is recognized in the device manager. (Advance remote check is required.)</p> <p>2) Reinstall the driver to the latest version. (Advance remote check is required.)</p> <p>3) Move the grabber to another PC slot and check if there is a problem.</p> <p>4) If the problem is not resolved, replace the grabber.</p> <p>5) If the same problem occurs after replacing the grabber, replace the PC.</p>
001	Appears when the frame grabber interface cannot be initialized, or memory for acquisition cannot be reserved	<p>1) Check whether the driver of the grabber is recognized in the device manager. (Advance remote check is required.)</p> <p>2) Reinstall the driver to the latest version. (Advance remote check is required.)</p> <p>3) Check the PC specification and the Windows &amp; Bios setting. (Advance remote check is required.)</p> <p>4) Move the grabber to another PC slot and check if the problem occurs.</p> <p>5) If the problem is not resolved, replace the grabber.</p> <p>6) If the same problem occurs after replacing the grabber, replace the PC.</p>

Tail code	Description	Solution
002	Appears when the MCU is not communicable, or the modem ring signal is in an improper state	1) Check whether the communication is available with the Main MCU. 2) If there is a problem, replace or reconnect the optic cable. 3) Replace the Main MCU.

### 11.2.3.2 Sub code – Capture Failure related error (020)

Tail code	Description	Solution
000	Appears when there is a capture error	1) Check the Windows & Bios setting. (Advance remote check is required.) 2) Check the grabber card. (Advance remote check is required.) ※ Check whether the device manager and grabber driver are recognizable. Check the FPGA version and conduct the pattern test. 3) Check whether the Dark is secured. (Advance remote check is required.) 4) If there is no problem, check whether 5V is applied to the sensor. 5) Connect the sensor and the PC Grabber with an optic cable to confirm whether the image data can be acquired. 5-1) If the equipment usually operates, check the optic cable and optic hub board for each section. ※ section: sensor~optic hub, optic hub~column, column~PC 6) If the problem is not resolved, replace the sensor.

### 11.2.3.3 Sub code – Reconstruction Failure related error (030)

Tail code	Description	Solution
001	Appears when bugs exist in	1) Check whether the graphics card is installed.

Tail code	Description	Solution
	VXM-file or there is insufficient memory	<p>2) Check whether the capturing PC satisfies the required PC specification (CPU, RAM, GPU).</p> <p>3) If the PC specification is usual, upgrade the graphics card driver to the latest version.</p> <p>4) If the problem is not resolved, replace the graphics card.</p>

**11.2.3.4 Subcode – Hardware related error (061)**

Tail code	Description	Solution
HW Error No	Appears when the error occurs during acquisition module operation	Follow the instructions stated in the hardware error code section according to the code number.

## 12. Cleaning and Maintenance

### NOTICE

The equipment must be installed and maintained on a flat surface.

### 12.1 Cleaning

### WARNING

Always turn off the power to the equipment and disconnect it from the power outlet before cleaning.

- Thoroughly clean the areas of the equipment that come in direct contact with the patient, such as the Chinrest and the Bite.
- Do not use spray cleaners or solvents as they could flow into the equipment and damage the electrical components or cause a fire.
- Do not use abrasive liquids such as acetone, gas, or oil, which may cause corrosion on the surface of the equipment.
- Do not use any cleaning products which contain silicon. They could potentially damage the equipment's electrical components.

The following table summarizes the standard cleaning procedures to be performed by the operator.

Components	Cleaning Process
Bite (Normal Bite, Special Bite A, and Special Bite B, Double Scan Support)	Clean with ethanol and gently wipe with a dry towel before the next patient.
Temple Supports	Clean with ethanol and gently wipe with a dry towel before the next patient.
Chinrest	Clean with ethanol and gently wipe with a dry towel before the next patient.
Computer and peripherals	Follow the manufacturers' instructions found in the accompanying manuals.
Outer covers of equipment	Wipe the unit with a dry cloth at the end of each day.

### IMPORTANT

Do not use cleaning agents in aerosol or spray form directly on the surface of the equipment.

## 12.2 Maintenance

**VATECH** requires periodic constancy tests to ensure image quality and safety for the patient and the operator.

Only **VATECH**-authorized technicians can perform inspection and service for the equipment. For technical assistance, contact the **VATECH** service center or your local **VATECH** representative.

### 12.2.1 Regular Maintenance



- Always turn off the equipment before performing any maintenance.
- Never remove equipment covers. There are no repairable parts inside.
- The only parts that can be replaced by the user are the input fuses, which must comply with the manufacturer's specifications.
- As a precaution against fire, the replacement should be one in the same type and range.

#### IMPORTANT

- There are no user-serviceable parts inside this equipment.
- If any service is required, please contact the **VATECH** service center or your local **VATECH** representative.

- Do not unplug cables by force.
- Do not expose the equipment or components in an area that is susceptible to water or humidity.
- Do not expose the equipment in an area of extreme fluctuation in temperature, poor ventilation, direct sunlight, dust, salt, etc.
- Keep all detachable components well organized and clean.
- Make sure that the equipment is well-grounded.
- Never try to modify this equipment, including the wires or cables. Doing so may damage it beyond repair.

### 12.2.2 Maintenance Task Checklist

Tasks	Period
Before the operation, ensure that the equipment is clean and ready for use. Make sure that all parts that come in direct contact with the patient have been cleaned thoroughly.	Daily
After using the equipment, make sure that the <b>Main Power Switch</b> has been turned off.	Daily
Ensure that the equipment is firmly plugged into a dedicated power source.	Daily
Ensure that the plug and the power cord are not heated abnormally.	Daily
Confirm that the LED indicator turns yellow when <b>the Exposure Switch</b> is pressed. Ensure that the LED indicator remains yellow for the entire duration of the exposure.	Daily
Ensure that the power cable is not kinked, broken, exposed and free of all other defects.	Daily
Confirm that activating the <b>Emergency Stop Switch</b> ceases the unit's operation. Pressing the <b>Emergency Stop Switch</b> should stop all movement of the equipment and X-ray emission.	Weekly
Ensure that all visible labels are intact and legible.	Weekly
Check for possible damages to the <b>Exposure Switch</b> cable.	Monthly
Confirm that the audio message is audible throughout the exposure.	Monthly

Left blank intentionally

## 13. Disposing of the Equipment

To reduce environmental contamination, this equipment is designed to be as safe as possible to use and dispose of. Many components of this equipment, except for some like the X-ray tube, are environment-friendly and can be recycled.

All parts and components which contain hazardous materials must be disposed of in accordance with disposal regulations (IEC 60601-1 6.8.2 j).

Parts	Materials	Recyclable	To the special disposal site	Hazardous waste; Needs Separate Collection
Frame and Covers	Aluminum and plastics	•		
Motors		•		
Circuit Boards		•		
Cables and Transformer	Copper	•		
	Steel	•		
	Oil		•	
Packing	Wood	•		
	Cardboard	•		
	Paper	•		
X-ray Tube				•
Sensor Head	Return the Sensor Head to <b>VATECH</b>			
Other parts			•	



This dental equipment shall not be disposed of as domestic garbage materials.

**IMPORTANT**

Clean / Disinfect / Sterilize the equipment before disassembling it and disposing of its parts.

**NOTICE**





Observe all regulations relevant to the disposal of waste in your country.

Left blank intentionally

# 14. Technical Specifications

## 14.1 Mechanical Specifications

### 14.1.1 Dimensions (unit = mm)



Item		Description
Weight	Without CEPH unit	134 kg (295.4 lbs. - without Base)
		187kg (412.3 lbs. – with Base)
	With CEPH unit	159 kg (350.5 lbs. - without Base)
		212 kg (467.4 lbs. - with Base)
Total Height	Without Base	Max. 2304 mm
	With Base	Max. 2335.5 mm
Dimensions during operation (Length x Width x Height)	Without CEPH unit	1488.7 mm (L) x 1125 mm (W) x 2304 mm (H) (without Base)
		1488.7 mm (L) x 1125 mm (W) x 2335.5 mm (H) (with Base)
	With CEPH unit	1488.7 mm (L) x 1874.1 mm (W) x 2304 mm (H) (without Base)
		1488.7 mm (L) x 1874.1 mm (W) x 2335.5 mm (H) (with Base)
Rotating Unit Vertical Movement		Max. 800 mm
Installation type		Base Stand / Wall Mount (Default: Wall Mount type)
Packing Box Organization		Main Box, CEPH Box (Optional), Base Box (Optional)

### 14.1.2 Image Magnification

Mode	FDD (mm)	FOD (mm)	ODD (mm)	Magnification
PANO	600	477.7	122.3	1 : 1.25
CEPH	1745	1524	221.0	1 : 1.14
CBCT	600	428.6	171.4	1 : 1.40

- **FDD:** Focal Spot to Detector Distance
- **FOD:** Focal Spot to Object Distance
- **ODD:** Object to Detector Distance (ODD = FDD - FOD)
- **Magnification** = FDD / FOD

## 14.2 Technical Specifications

### 14.2.1 X-ray Generator Specifications

#### Specifications


Item		Description	
Generator	Model	DG-07E22T2	
	Rated output power	1.6 kW	
	Inverter model name	INV-22	
	Type	Inverter	
	Normal/ Pulse	kVp	60 kV ~ 99 kV (1 kV increment)
		mA	4 mA ~ 16 mA (0.1 mA increment for CBCT, 1 mA increment for PANO and CEPH)
	Cooling	Cooling Protect (Optional fan cooling $\geq 60^{\circ}\text{C}$ )	
	Total filtration	Min. 2.5 mm Al	
	Default filtration	1.0 mm Al	
	Added filtration	1.5 mm Al (Fixed) / PANO and CEPH mode 1.5 mm Al (Fixed) + 3.0 mm Al (Automatically added) / CBCT mode	
Tube	Manufacturer	Canon Electron Tubes & Devices	
	Model	D-052SB (Stationary Anode type)	
	Focal spot size	0.5 mm (IEC 60336)	
	Target Angle	5 degree	
	Inherent Filtration	At least 0.8 mm Al equivalent at 50 kV	
	X-ray Coverage	95 mm x 380 mm at SID 550 mm	
	Anode Heat Content	35 kJ	
	Duty Cycle	1:60 or more (Exposure time : Interval time)	

**Test Condition**

Mode	Tube Voltage (kVp)	Tube Current (mA)	Exposure Time (s)
PANO	60 ~ 90	4 ~ 14	13.5
	60 ~ 90	4 ~ 14	11.5
	60 ~ 90	4 ~ 14	11.3
	60 ~ 90	4 ~ 14	11.1
	60 ~ 90	4 ~ 14	9.2
	60 ~ 90	4 ~ 14	7.7
	60 ~ 90	4 ~ 14	7.0
	60 ~ 90	4 ~ 14	6.8
	60 ~ 90	4 ~ 14	6.7
	60 ~ 90	4 ~ 14	6.1
	60 ~ 90	4 ~ 14	5.8
	60 ~ 90	4 ~ 14	5.7
	60 ~ 90	4 ~ 14	5.2
	60 ~ 90	4 ~ 14	5.0
	60 ~ 90	4 ~ 14	3.7
	60 ~ 90	4 ~ 14	3.5
	60 ~ 90	4 ~ 14	3.3
	60 ~ 90	4 ~ 14	2.8
	CEPH	60 ~ 99	4 ~ 16
60 ~ 99		4 ~ 15	2.4
60 ~ 99		4 ~ 15	3.9
60 ~ 99		4 ~ 14	4.9
60 ~ 99		4 ~ 14	5.4
CBCT	60 ~ 99	4 ~ 12	9.0
	60 ~ 99	4 ~ 12	4.9


### Maximum Rating Charts

Constant potential high-voltage generator  
Nominal Focus Spot Value: 0.5




### Emission & Filament Characteristics


Constant potential high-voltage generator  
Nominal Focus Spot Value: 0.5



### ■ Anode Thermal Characteristics



### ■ X-ray Housing Assembly Tube Characteristics





### 14.2.2 Detector Specifications

Item	Description	
	PANO & CBCT	CEPH
Model	Xmaru1314CF	Xmaru2602CF
Detector Type	CMOS photodiode array	
Pixel size	198 $\mu\text{m}$ @ 4X4 binning (49.5 $\mu\text{m}$ @ no binning)	200 $\mu\text{m}$ @ 2x2 binning (100 $\mu\text{m}$ @ no binning)
Active area	CBCT - 127.5 mm X 135.8 mm PANO - 5.9 mm X 135.8 mm	15.6 mm x 259 mm
Frame Rate	~ 108 Hz @ 4X4 binning	~ 330 Hz @ 2x2 binning
Analogue-Digital Conversion	14 bits	
Operating condition	10 ~ 35 $^{\circ}\text{C}$ (Temperature) 10 ~ 75 % (Humidity)	
Storage condition	-10 ~ 60 $^{\circ}\text{C}$ (Temperature) 10 ~ 75 % (Humidity)	
Sensor size	165 mm (L) x 230 mm (W) x 27 mm (H)	110 mm (L) x 279 mm (W) x 20 mm (H)
Sensor weight	1450 g	1050 g
Converter	CsI : Ti	
Energy Range	50 ~ 120 kVp	
Readout	Charge amplifier array	
Video Output	Optic	
MTF	> 45 % @ 1 lp/mm, > 10 % @ 2.5 lp/mm	> 40 % @ 1 lp/mm > 8 % @ 2.5 lp/mm
DQE	> 60 % @ ~0 lp/mm	> 70 % @ ~0 lp/mm
Dynamic Range	> 80 dB	> 70 dB

### 14.3 Electrical Specifications

Item	Description
Power supply voltage	100 - 240 V ~
Frequency	50 / 60 Hz
Power rating	2.0 kVA
Accuracy	Tube Voltage (kVp) $\pm$ 10 %, Tube Current (mA) $\pm$ 20 %, Exposure Time (s) $\pm$ (5 % + 50 ms)

- The input line voltage depends on the local electrical distribution system.
- Allowable input voltage fluctuation requirement:  $\pm$ 10 %.
- Mode of operation: Continuous operation with intermittent loading - Needs waiting time (at least 60 times the exposure time) before the next exposure begins.
- Column operation time: Max. 2 min. On / 18 min. Off (Ratio 1:9)



**NOTICE**

- To assure line voltage quality, a separate 3-core grounded power cable connected directly to the central distribution panel with an over-current circuit breaker rated for 20A must be used.
- Maximally allowed deviation of the tube voltage / tube current / exposure time:  
Tube Voltage (kVp)  $\pm 10\%$  / Tube Current (mA)  $\pm 20\%$  / Exposure Time (s)  $\pm (5\% + 50\text{ ms})$  according to IEC 60601-2-63.
- The mains resistance should not exceed 0.045 ohms at 100 V and 0.19 ohm at 240 V.

**14.4 Environmental Specifications**

	Item	Description
During Operation	Temperature	10 ~ 35 °C
	Relative humidity	30 ~ 75 %
	Atmospheric pressure	860 ~ 1060 hPa
During Transport and Storage	Temperature	-10 ~ 60 °C
	Relative humidity	10 ~ 75 %
	Atmospheric pressure	860 ~ 1060 hPa

Left blank intentionally

## 15. Appendices

### 15.1 Recommended X-ray Exposure Tables

#### 15.1.1 PANO Mode

##### Exposure Condition

Mode	Image Option	Gender / Age group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
PANO Examination	High Resolution	Man	Hard	75	12
			Normal	74	12
			Soft	73	12
		Woman	Hard	74	12
			Normal	73	12
			Soft	72	12
		Child	Hard	68	10
			Normal	67	10
			Soft	66	10
	Green (Optional)	Man	Hard	75	14
			Normal	74	14
			Soft	73	14
		Woman	Hard	74	14
			Normal	73	14
			Soft	72	14
Child		Hard	68	12	
		Normal	67	12	
		Soft	66	12	
PANO Examination (Insight PAN) (Optional)	High Resolution	Man & Woman	Normal	70	7
		Child	Normal	67	7
SPECIAL Examination	N/A	Man	Hard	75	14
			Normal	74	14
			Soft	73	14
		Woman	Hard	74	14
			Normal	73	14
			Soft	72	14
		Child	Hard	68	12
			Normal	67	12
			Soft	66	12

**Scan Time / Exposure Time**

Examination Mode	Arch Type	Examination Mode	High Resolution		Green (Optional)		
			Scan Time (s)	Exposure Time (s)	Scan Time (s)	Exposure Time (s)	
PANO Examination	Narrow	Standard	14.1	13.5	7.0	7.0	
		Right	14.1	6.8	7.0	3.5	
		Front	14.1	11.3	7.0	5.8	
		Left	14.1	6.8	7.0	3.5	
	Normal	Standard	14.1	13.5	7.0	7.0	
		Right	14.1	6.8	7.0	3.5	
		Front	14.1	11.3	7.0	5.8	
		Left	14.1	6.8	7.0	3.5	
	Wide	Standard	14.1	13.5	7.0	7.0	
		Right	14.1	6.8	7.0	3.5	
		Front	14.1	11.3	7.0	5.8	
		Left	14.1	6.8	7.0	3.5	
	Child	Standard	12.0	11.5	6.8	6.7	
		Right	12.0	5.7	6.8	3.3	
		Front	12.0	9.2	6.8	5.2	
		Left	12.0	5.7	6.8	3.3	
	Orthogonal	Standard	14.1	13.5	7.0	7.0	
		Right	14.1	6.7	7.0	3.5	
		Front	14.1	11.1	7.0	5.7	
		Left	14.1	6.7	7.0	3.5	
		Bitewing	14.1	9.2	7.0	5.0	
		Bitewing Incisor (Optional)	14.1	2.8	7.0	1.4	
		Bitewing Right	14.1	5.0	7.0	2.8	
		Bitewing Left	14.1	5.0	7.0	2.8	
	SPECIAL Examination	-	TMJ LAT Open	14.1	6.7	14.1	6.7
			TMJ LAT Close				

Examination Mode	Arch Type	Examination Mode	High Resolution		Green (Optional)	
			Scan Time (s)	Exposure Time (s)	Scan Time (s)	Exposure Time (s)
		TMJ PA Open (Optional)	10.0	6.1	10.0	6.1
		TMJ PA Close (Optional)				
		Sinus LAT (Optional)	4.0	3.7	4.0	3.7
		Sinus PA	8.8	7.7	8.8	7.7

- Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.
- Exposure Time: The actual time that the patient is exposed to the X-ray emission.

## 15.1.2 CEPH Mode

Exposure Condition

Examination Program	Image Option	Gender / Age group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
Lateral	High Resolution	Man	Hard	92	15.0
			Normal	90	15.0
			Soft	88	15.0
		Woman	Hard	90	15.0
			Normal	88	15.0
			Soft	86	15.0
		Child	Hard	88	15.0
			Normal	86	15.0
			Soft	84	15.0
	Green	Man	Hard	92	16.0
			Normal	90	16.0
			Soft	88	16.0
		Woman	Hard	90	16.0
			Normal	88	16.0
			Soft	86	16.0
Child		Hard	88	16.0	
		Normal	86	16.0	
		Soft	84	16.0	
Full Lateral (Optional)	High Resolution / Green	Man	Hard	92	14.0
			Normal	90	14.0
			Soft	88	14.0
		Woman	Hard	90	14.0
			Normal	88	14.0
			Soft	86	14.0
		Child	Hard	88	14.0
			Normal	86	14.0
			Soft	84	14.0

Examination Program	Image Option	Gender / Age group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
PA SMV Waters' view	High Resolution	Man	Hard	92	14.0
			Normal	90	14.0
			Soft	88	14.0
		Woman	Hard	90	14.0
			Normal	88	14.0
			Soft	86	14.0
		Child	Hard	88	14.0
			Normal	86	14.0
			Soft	84	14.0
	Green	Man	Hard	92	15.0
			Normal	90	15.0
			Soft	88	15.0
		Woman	Hard	90	15.0
			Normal	88	15.0
			Soft	86	15.0
Child		Hard	88	15.0	
		Normal	86	15.0	
		Soft	84	15.0	
Carpus	High Resolution / Green	Man	Hard	90	6.0
			Normal	88	6.0
			Soft	86	6.0
		Woman	Hard	88	6.0
			Normal	86	6.0
			Soft	84	6.0
		Child	Hard	86	6.0
			Normal	84	6.0
			Soft	82	6.0

**Scan Time / Exposure Time**

Examination Program	High Resolution		Green	
	Scan Time (s)	Exposure Time (s)	Scan Time (s)	Exposure Time (s)
Lateral	3.9	3.9	1.9	1.9
Full Lateral (Optional)	5.4	5.4	3.9	3.9
PA	4.9	4.9	2.4	2.4
SMV	4.9	4.9	2.4	2.4
Waters' view	4.9	4.9	2.4	2.4
Carpus	4.9	4.9	2.4	2.4

- Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.
- Exposure Time: The actual time that the patient is exposed to the X-ray emission.

### 15.1.3 CBCT Mode

#### Exposure Area

FOV (cm)	Vertical Position	Horizontal Position		
		Right	Center	Left
16x9	Occlusion	X	O	X
12x9	Occlusion	X	O	X
	TMJ	O	X	O
	Airway	X	O	X
8x9	Occlusion	O	O	O
	TMJ	O	X	O
8x5	Maxilla / Mandible	O	O	O
5x5	Maxilla / Mandible	Right Molar / Right / Incisor / Left / Left Molar		
Double scan (Optional)	First Exposure: Mandible Second Exposure: Maxilla	X		

#### Exposure Condition

FOV (cm)	Image Option	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
16x9 & 12x9 & Double scan (Optional)	High Resolution	Man	Hard	95	8.0
			Normal	94	8.0
			Soft	93	8.0
		Woman	Hard	95	7.7
			Normal	94	7.7
			Soft	93	7.7
		Child	Hard	95	7.4
			Normal	94	7.4
			Soft	93	7.4
	Green	Man	Hard	88	6.1
			Normal	87	6.1
			Soft	86	6.1
Woman		Hard	88	5.8	
		Normal	87	5.8	

FOV (cm)	Image Option	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)	
8x9, 8x5, 5x5.		Child	Soft	86	5.8	
			Hard	88	5.5	
			Normal	87	5.5	
				Soft	86	5.5
	High Resolution	Man	Hard	95	8.0	
			Normal	94	8.0	
			Soft	93	8.0	
		Woman	Hard	95	7.7	
			Normal	94	7.7	
			Soft	93	7.7	
		Child	Hard	95	7.4	
			Normal	94	7.4	
			Soft	93	7.4	
		Green	Man	Hard	88	6.1
				Normal	87	6.1
Soft	86			6.1		
Woman	Hard		88	5.8		
	Normal		87	5.8		
	Soft		86	5.8		
Child	Hard		88	5.5		
	Normal		87	5.5		
	Soft		86	5.5		

**Scan Time / Exposure Time**

FOV (cm)	Scan Time (s) (High Resolution / Green)	Exposure Time (s) (High Resolution / Green)
16x9	9.0	9.0
12x9	9.0	9.0
8x9	4.9	4.9
8x5	4.9	4.9
5x5	4.9	4.9
Double scan (Optional)	18.0	18.0

- Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.
- Exposure Time: The actual time that the patient is exposed to the X-ray emission.

**Reconstruction Time / File Size (Measured Object: Skull)**

FOV (cm)	Voxel Size (mm)	Reconstruction Time (s)		File Size (MB)
		Fast	HD	
Double scan (Optional)	0.3	60	-	282
16x9	0.2	22	145	550
	0.3	15	78	163
12x9	0.2	16	101	309
	0.3	13	58	92
8x9	0.12	18	160	635
	0.2	11	56	138
8x5	0.12	12	98	359
	0.2	9	37	77
5x5	0.08	18	114	466
	0.12	8	44	138

- The above data is obtained from a computer system which is based on Intel E5-1607 v3@3.10GHz (16GB of RAM) and NVIDIA GeForce GTX1060 6GB.
- Image reconstruction time varies depending on computer specifications and working conditions.

### 15.1.4 3D MODEL Scan Mode

#### Exposure Area

FOV (cm)	MODEL Type	Horizontal Position		
		Right	Center	Left
8x9	Upper (Maxilla)	X	O	X
	Lower (Mandible)	X	O	X

#### Exposure Condition

FOV (cm)	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
8x9	Man / Woman / Child	Hard / Normal / Soft	95	8.0

#### Scan Time / Exposure Time

FOV (cm)	Scan Time (s)	Exposure Time (s)
8x9	9.0	9.0

- Scan Time: The actual time that the equipment shoots the plaster cast except for the initial acceleration and late deceleration stages.
- Exposure Time: The actual time that the plaster cast is exposed to the X-ray emission.

#### Reconstruction Time / File Size

FOV (cm)	Voxel Size (mm)	Reconstruction Time (s)	File Size (MB)
8x9	0.12	269	635

- The above data is obtained from a computer system which is based on Intel i7-6700 and NVIDIA GeForce GTX1060 6GB.
- Image reconstruction time varies depending on computer specifications and working conditions.

## 15.2 X-ray Dose Data

### 15.2.1 DAP (Dose Area Product)

The X-ray dose data is extracted from the X-ray Dose Test Report for **Green16 (PHT-65LHS)**.

X-ray Dose Test Report for the **Green16 (PHT-65LHS)** maintains dosimetric evaluation that the **VATECH** dental diagnostic system meets all requirements specified in the IEC Collateral Standard. To limit unnecessary exposure to the patient, operator or other staff, **Green16 (PHT-65LHS)** is designed to comply with IEC 60601-1-3 Part 1 General Requirements for Safety.

Test Hardware	
Brand Name (Model)	<b>Green16 (PHT-65LHS)</b>
Sensor Type	PANO & CBCT: Xmaru1314CF CEPH: Xmaru2602CF
X-ray Generator	DG-07E22T2
Tube	D-052SB

DAP (Dose Area Product) is a quantity used in assessing the radiation risk from diagnostic X-ray examination procedures. It is defined as the absorbed dose multiplied by the area irradiated, expressed in gray square centimeters ( $\text{mGy}\cdot\text{cm}^2$ ). Despite the limitation, DAP is the best way to predict effective dose value and currently the most convenient method for patient doses monitoring.


- *PHE (Public Health England) recommends that any national reference dose that achievable dose (DAP) value of 250 [ $\text{mGy}\cdot\text{cm}^2$ ] for a clinical protocol for a standard male patient.*

**DAP (Dose Area Product) Calculation**

$$\text{DAP}[\text{mGy} \cdot \text{cm}^2] = \text{Dose}[\text{mGy}] \times \text{Exposed Area}[\text{cm}^2]$$

**NOTICE**

When you need more information on DAP measurement procedures or test results for the equipment, please contact **VATECH** service center or your local **VATECH** representative and get assistance from **VATECH**-authorized technicians

**Measurement Overview****Results**

Mode	Exposure Condition	DAP [mGy·cm <sup>2</sup> ]
PANO Adult Man Normal (High Resolution)	74 kVp / 12.0 mA / 13.5 s	110
PANO Child Normal (High Resolution)	67 kVp / 10.0 mA / 11.5 s	55
PANO Adult Man Normal (Green)	74 kVp / 14.0 mA / 7.0 s	65
PANO Child Normal (Green)	67 kVp / 12.0 mA / 6.7 s	35
CEPH Adult Man LAT (High Resolution)	90 kVp / 15.0 mA / 3.9 s	41
CEPH Child LAT (High Resolution)	86 kVp / 10.0 mA / 3.9 s	38
CEPH Adult Man LAT (Green)	90 kVp / 16.0 mA / 1.9 s	25
CEPH Child LAT (Green)	86 kVp / 16.0 mA / 1.9 s	23
CBCT 16x9 Adult Man (High Resolution)	94 kVp / 8.0 mA / 9.0 s	1435
CBCT 16x9 Adult Man (Green)	87 kVp / 6.1 mA / 9.0 s	913
CBCT 12x9 Adult Man (High Resolution)	94 kVp / 8.0 mA / 9.0 s	1147
CBCT 12x9 Adult Man (Green)	87 kVp / 6.1 mA / 9.0 s	730
CBCT 8x9 Adult Man (High Resolution)	94 kVp / 8.0 mA / 4.9 s	802

Mode	Exposure Condition	DAP [mGy·cm <sup>2</sup> ]
CBCT 8x9 Adult Man (Green)	87 kVp / 6.1 mA / 4.9 s	510
CBCT 8x5 Maxilla/Mandible (High Resolution)	94 kVp / 8.0 mA / 4.9 s	490
CBCT 8x5 Maxilla/Mandible (Green)	94 kVp / 8.0 mA / 4.9 s	307
CBCT 5x5 Maxilla Adult Man (High Resolution)	94 kVp / 8.0 mA / 4.9 s	331
CBCT 5x5 Maxilla Adult Man (Green)	87 kVp / 6.1 mA / 4.9 s	211
CBCT 5x5 Mandible Adult Man (High Resolution)	94 kVp / 8.0 mA / 4.9 s	331
CBCT 5x5 Mandible Adult Man (Green)	87 kVp / 6.1 mA / 4.9 s	211
CBCT Double Scan Adult Man (High Resolution) (Optional)	94 kVp / 8.0 mA / 9.0 s	2870
CBCT Double Scan Adult Man (Green) (Optional)	87 kVp / 6.1 mA / 9.0 s	1826

- In Double Scan mode, exposure the FOV 16x9 twice, then stitch the acquired images to provide one image.


## 15.2.2 Leakage Dose

X-ray Leakage Dose Test is performed to protect patients against excessive and unnecessary radiation that is not purposed, and this document evaluates leakage dose amount based on the following standard defined by IEC regulation and has been performed by covering each collimator region in use.

### Standard

National Deviation	Terminology	Permissive Range
International Standard IEC 60601-1-3	Leakage	limits leakage at 1M from the source to 100 mR in 1hr

### Measurement Overview




### Test Condition

Test mode	Tube Voltage (kVp)	Tube Current (mA)	Exposure Time (s)
PANO Adult Man Normal (High Resolution)	99	14.0	13.5
CEPH Adult Man LAT (High Resolution)	99	16.0	3.9
CBCT16x9 Adult Man (High Resolution)	99	12.0	9.0

- In Double Scan mode, exposure the FOV 16x9 twice, then stitch the acquired images to provide one image.

**Results**


Direction [°]	PANO Adult Man Normal (High Resolution) [mR/h]	CEPH Adult Man LAT (High Resolution) [mR/h]	CBCT 16x9 Adult Man (High Resolution) [mR/h]
0	7.3	5.4	14.6
45	5.8	4.5	8.3
90	1.1	2.8	4.8
135	3.5	4.5	5.2
180	10.5	21.0	11.7
225	22.4	11.8	33.5
270	6.8	18.5	7.4
315	24.5	5.5	48.7



### 15.2.3 Scattered Dose

X-ray Scattered Dose data concerning different angle and distance is examined for recommendations about appropriate radiation level insignificant zones of occupancy and the effectiveness of protective shielding facility around the patient's position. This information states the identity and intended position of the tested phantom and scattered dosimetric evaluation under the defined scope and test circumstances to ensure the magnitude of risks to the operator and staff, during both accident situations and routine work.

#### Measurement Overview



#### Test Condition


Test mode	Tube Voltage (kVp)	Tube Current (mA)	Exposure Time (s)
PANO Adult Man Normal (High Resolution)	99	14.0	13.5
CBCT 16x9 Adult Man (High Resolution)	99	12.0	9.0
CEPH Adult Man LAT (High Resolution)	99	16.0	3.9

- In Double Scan mode, exposure the FOV 16x9 twice, then stitch the acquired images to provide one image.

## Results

### ■ (PANO / Adult, Man / Normal / High Resolution)

Direction [°]	1 m [mR]	1.5 m [mR]	2 m [mR]
0	0.079	0.049	0.021
45	0.039	0.025	0.016
90	0.079	0.050	0.020
135	0.050	0.032	0.013
180	0.162	0.085	0.034
225	0.175	0.082	0.035
270	0.272	0.077	0.031
315	0.224	0.066	0.030



■ CBCT / 16x9 / Adult, Man / High Resolution


Direction [°]	1 m [mR]	1.5 m [mR]	2 m [mR]
0	1.727	0.971	0.482
45	1.817	0.952	0.499
90	0.656	0.170	0.044
135	2.392	1.121	0.525
180	2.541	1.104	0.480
225	2.399	1.095	0.500
270	2.415	1.185	0.581
315	1.874	0.936	0.467



### ■ CEPH / Full Lateral / Adult, Man / High Resolution

Direction [°]	1 m [mR]	1.5 m [mR]	2 m [mR]
0	0.031	0.021	0.012
45	0.013	0.009	0.003
90	0.011	0.007	0.004
135	0.033	0.006	0.004
180	0.025	0.013	0.004
225	0.045	0.034	0.019
270	-	-	0.022
315	0.073	0.048	0.025

- Data of 1 m and 1.5 m at 270 ° are not measured since the Ion chamber is located between the generator and the object.



## 15.3 Electromagnetic Compatibility (EMC) Information

### Guidance and manufacturer's declaration - electromagnetic emissions

The **PHT-65LHS** is intended for use in the electromagnetic environment specified as below. The customer or the user of the **PHT-65LHS** should assure that it is used in such an environment.

Immunity test	Compliance	Electromagnetic environment - Guidance
RF emissions CISPR 11	Group 1	The <b>PHT-65LHS</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The <b>PHT-65LHS</b> is suitable for use in all establishments other than domestic and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	<p><b>Warning:</b> This equipment/system is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the <b>PHT-65LHS</b> or shielding the location.</p>

(NOTE) It is essential that the actual RF shielding effectiveness and filter attenuation of the shielded location be verified to ensure that they meet or exceed the specified minimum values.

**Guidance and manufacturer's declaration - electromagnetic immunity**

The **PHT-65LHS** is intended for use in the electromagnetic environment specified below. The customer or the user of the **PHT-65LHS** should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV Contact  ± 2 kV, ± 4 kV, ± 8 kV, ±15 kV air	± 8 kV Contact  ± 2 kV, ± 4 kV, ± 8 kV, ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	The main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U <sub>T</sub> : 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U <sub>T</sub> : 1 cycle and 70 % U <sub>T</sub> : 25/30 cycles Single-phase: at 0°	0 % U <sub>T</sub> : 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U <sub>T</sub> : 1 cycle and 70 % U <sub>T</sub> : 25/30 cycles Single-phase: at 0°	The main power quality should be that of a typical commercial or hospital environment. If the user of the <b>PHT-65LHS</b> image intensifier requires continued operation during main power interruptions, it is recommended that the <b>PHT-65LHS</b> image intensifier be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U<sub>T</sub> is the a.c. mains voltage before application of the test level.

**Guidance and manufacturer’s declaration - electromagnetic immunity**

The **PHT-65LHS** is intended for use in the electromagnetic environment specified below. The customer or the user of the **PHT-65LHS** should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	V <sub>1</sub> =3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the <b>PHT-65LHS</b> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	E <sub>1</sub> =3V/m	<p><b>Recommended separation distance:</b></p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup></p> <p>Interference may occur in the vicinity of the equipment marked with the following symbol:</p> 

NOTE 1) At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

The **PHT-65LHS** is intended for use in the electromagnetic environment specified below. The customer or the user of the **PHT-65LHS** should assure that it is used in such an environment.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location where the **PHT-65LHS** is used exceeds the applicable RF compliance level above, the **PHT-65LHS** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **PHT-65LHS**.

<sup>b</sup> Over the frequency range 150kHz to 80MHz, field strengths should be less than [V<sub>1</sub>] V/m.

### **Recommended separation distances between portable and mobile RF communications equipment and the PHT-65LHS**

This is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **PHT-65LHS** can help Prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **PHT-65LHS** as recommended below, according to the maximum output power of the communications equipment.

#### **Separation distance according to the frequency of transmitter [m]**

##### **IEC 60601-1-2**

Frequency of Transmitter	150kHz to 80MHz	80MHz to 800MHz	800MHz to 2.7GHz
Equation	$d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$	$d = \left[ \frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
The rated maximum output power of the transmitter [W]	V <sub>1</sub> =3Vrms Separation Distance (meters)	E <sub>1</sub> =3V/m Separation Distance (meters)	E <sub>1</sub> =3V/m Separation Distance (meters)
0.01	0.116	0.1166	0.2333
0.1	0.368	0.3687	0.7378
1	1.166	1.1660	2.3333
10	3.687	3.6872	7.3785
100	11.660	11.6600	23.333

This is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **PHT-65LHS** can help Prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **PHT-65LHS** as recommended below, according to the maximum output power of the communications equipment.

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a. Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **PHT-65LHS** is used exceeds the applicable RF compliance level above, the **PHT-65LHS** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **PHT-65LHS**
- b. Over the frequency range 150kHz to 80MHz, field strengths should be less than  $[V_1]$  V/m.

### Immunity and Compliance Level

Immunity test	IEC 60601-1-2 Test level	Actual Immunity Level	Compliance Level
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	3Vrms
Radiated RF IEC 61000-4-3	3Vrms 80MHz to 2.7GHz	3V/m	3V/m


## 15.4 Hand-wrist Image Evaluation References

**Radiographic Evaluation of Skeletal Maturation. A clinically oriented method based on hand-wrist films.**

**Fishman LS. 1982**

### **The system of Skeletal Maturation Assessment (SMA)**

The System uses only four stages of bone maturation; all found at six anatomical sites located on the thumb, third finger, fifth finger, and radius, as seen in Fig.1. Eleven discrete adolescent skeletal maturational indicators (SMI's), covering the entire period of adolescent development, are found on these six sites (Fig.1 orange circles).



[Fig1. The site of skeletal maturity indicators]

## Skeletal Maturity Indicators (SMI)

A system of skeletal maturation assessment based on four stages of bone maturation at six anatomical sites in the hand wrist.




Fig. 2 Radiographic identification of skeletal maturity indicators.  
 A. Epiphysis equal in width to diaphysis.  
 B. Appearance of adductor sesamoid of the thumb.  
 C. Capping of epiphysis.  
 D. Fusion of epiphysis.

[Fig2. Radiographic identification of skeletal maturity indicators]

### A. The width of epiphysis as wide as the diaphysis

1. Third finger – a Proximal phalanx
2. Third finger – a middle phalanx
3. Fifth finger – a middle phalanx

### B. Ossification

1. Adductor sesamoid of thumb

### C. Capping of epiphysis


1. Third finger – a distal phalanx
2. Third finger – a middle phalanx
3. Fifth finger – a middle phalanx

#### D. Fusion

1. Third finger – a distal phalanx
2. Third finger – a proximal phalanx
3. Third finger – a middle phalanx
4. Radius

#### Eleven Skeletal maturity indicators (SMIs)


The System uses only four stages of bone maturation; all found at six anatomical sites located on the thumb, third finger, fifth finger, and radius, as seen in Fig.1. Eleven discrete adolescent skeletal maturational indicators (SMI's), covering the entire period of adolescent development, are found on these six sites (Fig.1 orange circles).



[Fig3. Eleven Skeletal maturity indicators (SMIs)]

A systematic observational scheme such as that shown in the figure above can further facilitate SMI evaluation. With this approach, key stages are checked first, rather than looking for maturity indicators in numerical order. A useful step is to determine if the adductor sesamoid of the thumb can be seen (orange circle). If not, then the applicable SMI will be one of those associated with early epiphyseal widening rather than capping. If the sesamoid is visible, then either the sesamoid or an SMI based on capping or fusion will be applicable.

## HAND-WRIST OBSERVATION SCHEME



## 15.5 Acquiring Images for Pediatric Dental Patients

### 15.5.1 Age Group: Classification Table

Ages are classified loosely into the following correspondence between FDA definition and one used in this manual.

Age Group	FDA's standard	VATECH's Standard
Infant	1 month to 2 years	N/A
Child	2 ~ 12 years of age	Child
Adolescent	12 ~16 years of age	Adult
Other	16 ~ 21 years of age	
Adult	> 21 years of age	

### 15.5.2 Positioning the Pediatric Dental Patients

1. Use a laser light beam guide to locate the midsagittal plane. Direct patient focusses on mirroring reflection. Affix decal to mirror to aid the patient in maintaining the correct position throughout the exposure.
2. Move the Chinrest into a position that is slightly higher than the patient's chin height before requesting that the weak place chin onto the rest. Direct the patient to assume a position that resembles the erect stance of a soldier.
3. Direct the patient to stick out the chest while dropping the chin down. While holding the unit handles for stability, direct the patient to take a half step in toward the vertical column of the X-ray device into a position that feels as if he/she is slightly leaning backward.
4. Direct the patient to close lips around the Bite Block during the exposure.
5. Direct the patient to swallow and note the flat position of the tongue. Request that the patient sucks in the cheeks, pushing the tongue into the correct flat position against the palate and maintain this position throughout the exposure.

### <How to product error-free radiographic images for the pediatric patient>

(<http://www.dimensionsofdentalhygiene.com/print.aspx?id=3612>)

- By Evelyn M. Thomson, BSDH, MS

Panoramic radiographs are often recommended for assessing growth and development of the pediatric patient and for evaluation of developing third molars during adolescence.<sup>1-3</sup> While the panoramic technique seems relatively straightforward, producing a diagnostic quality image of the pediatric patient requires a mastery of technical skill.<sup>4</sup> Modern panoramic x-ray equipment is designed for ease of use, yet studies continue to demonstrate a high incidence of errors.<sup>5-7</sup> Positioning errors may occur at an even higher rate in pediatric panoramic radiographs.<sup>7</sup> The goal of the dental hygienist is to maximize the use of panoramic imagery in the assessment of the pediatric patient while minimizing the occurrence of retakes that result from the radiographic error.

#### Producing A Quality Panoramic Image

A quality panoramic radiograph should image all of the teeth, erupted and unerupted, in both the maxillary and mandibular arches from condyle to condyle in the horizontal dimension, and from the superior third of the orbit in the superior region to the inferior border of the mandible in the inferior region.<sup>8,9</sup> The arches should appear straight or slightly U-shaped with the occlusal plane parallel to the horizontal edges of the film **(Figure 1)**. The anterior teeth must not be magnified or diminished in size and overlapping of adjacent posterior teeth should be kept to a minimum.



**Figure 1:** Example of a diagnostically acceptable panoramic radiograph of an adolescent patient undergoing orthodontic intervention. (Courtesy of Jamie Mace and Will Wright of Schick Technologies Inc.)

The most critical component in producing a diagnostically acceptable panoramic image is patient positioning. All panoramic x-ray machines have guidelines to assist with

positioning the dental arches within the three dimensions of the focal trough, an area where the anatomical structures will be imaged in relative clarity. Most panoramic x-ray machines have a bite block to indicate the correct anterior-posterior position, or how far forward or back the patient should be positioned, side positioner guides for determining the correct lateral alignment, and chin rest to correctly locate the superior-inferior dimension or how far up or down the chin should be positioned.<sup>4,10</sup> Panoramic x-ray machines are available with a mirror and laser light beam guide that shines on the patient's face to illustrate various anatomical planes (**Figure 2**). Incorrectly positioning the patient in any of these three dimensions will produce unique and distinct radiographic image errors (**Table 1**).



**Figure 2:** Laser light beam guides that assist with determining correct patient positioning.

**Table 1.** Common Panoramic Positioning Errors

Error	Cause	Corrective action	Tips for pediatric patients
Anterior teeth narrow Severe posterior overlap Vertebrae superimposed over condyles	Arches positioned too far anterior	Position anterior teeth in appropriate position on bite guide.	Use a cotton roll to fill in missing primary teeth or partially erupted permanent teeth. Adapt adult recommendation for direction of laser light beam guide for use with primary teeth.
Anterior teeth wide, blurred out of image Condyles not imaged	Arches positioned too far posterior	Locate appropriate position with anterior laser light guide.	Observe laser light beam guide on both the right and left sides.
Teeth on the right side appear narrowed, severely overlapped Teeth on the left side appear broad, poorly defined Condyles asymmetrical in width and height	Arches tipped or tilted to the right	Position the midsagittal plane perpendicular to the floor.	Use laser light beam guide to locate midsagittal plane. Direct patient focus to mirror reflection. Affix decal to mirror to aid patient in maintaining the correct position throughout exposure.
Teeth on the left side appear narrowed, severely overlapped Teeth on the right side appear broad and poorly defined Condyles asymmetrical in width and height	Arches tipped or tilted to the left		
Flat, downward-turned, "frown" appearance to the occlusal plane Palate appears as a widened, thick, dense radiopacity Condyles flare out off the edges of the image Anterior teeth appear wide, elongated	Arches positioned too far superior	Position the Frankfort or the canthomeatal plane parallel to the floor, or the ala-tragus line 5° down toward the floor.	Move chin rest into a position that is slightly higher than the patient's chin height before requesting that the patient place chin onto the rest. Direct the patient to assume a position that resembles the erect stance of a soldier.
Exaggerated upward curve of the occlusal plane creating a "smile" appearance Hyoid bone superimposed over the mandible Condyles tilt inward Anterior teeth appear narrowed; elongated in the maxilla and foreshortened in the mandible	Arches positioned too far inferior		
Pyramid-shaped radiopacity superimposed over the anterior teeth	Patient in slumped position	Position the back and neck straight.	Direct the patient to stick out the chest while dropping the chin down. While holding the unit handles for stability, direct the patient to take a half step in toward the vertical column of the x-ray machine into a position that feels as if he/she is slightly leaning backward.
Radiolucent shadow of the commissure superimposed over the teeth, mimicking caries	Lips not closed around bite block	Position the lips around the bite block.	Direct the patient to keep the lips closed around the bite block during the exposure.
Radiolucency superimposed over the maxillary teeth apices	Tongue not placed against palate	Position the tongue flat against the roof of the mouth.	Direct the patient to swallow and note the flat position of the tongue. Request that the patient suck in the cheeks, pushing the tongue into the correct flat position against the palate and maintain this position throughout the exposure.

FEBRUARY 2009

Dimensions OF DENTAL HYGIENE

### Anterior-Posterior Positioning Error

When the arches are positioned incorrectly in the anterior-posterior direction, distortion or ghosting of the anterior anatomy occurs. Unerupted teeth in the anterior region may not be imaged on the radiograph if positioned outside of the focal trough. It is important to note that an error of only 3 mm to 4 mm in either direction will result in a significantly compromised image.<sup>11</sup> When the arches are positioned too far anterior, the anterior teeth will appear narrow and diminished in size. The vertebrae of the spinal column may be superimposed over the condyles at the edges of the film and, depending on the size of the child, may be superimposed over the rami of the mandible blocking a clear view of the posterior teeth (**Figure 3**). When the arches are positioned too far posteriorly, the anterior teeth will appear broad or widened. If the position is excessively posterior, anterior teeth may be blurred entirely from the image and the condyles may be cut off from the edges of the film.



**Figure 3:** Incorrect position too far anteriorly. Note the narrow anterior teeth and superimposition of the spinal column over the condyles. The radiolucency superior to the maxillary apices indicates that the tongue was not placed against the palate. An open lip line can also be detected.

To avoid these imaging errors, the anterior teeth must occlude edge-to-edge onto the designated area of the bite block. Achieving this position is easily compromised during exfoliation of primary teeth, making precise occlusion difficult when one tooth or multiple teeth are missing or partially erupted. A cotton roll may be attached to the bite block to fill in the space created by the missing tooth or teeth. Additionally, an adjustment may be necessary when using a laser light beam guide. The manufacturer's instructions for directing the laser light beam at a predetermined tooth or interproximal space usually apply to adult patients. These instructions may need to be modified for the pediatric patient with primary or mixed dentition.

### Lateral Left-Right Positioning Error

When the arches are positioned incorrectly in the lateral left-right dimension, the posterior teeth on one side will appear broad or widened, while the teeth on the other side will appear narrowed or diminished in width and severely overlapped (**Figure 4**). This image distortion is like that which occurs with an incorrect anterior-posterior position. When the arches are rotated or tilted, the posterior teeth on one side move out of the focal trough to a position further away (back) from the image receptor, while the opposite side simultaneously moves closer (forward) to the image receptor. Depending on the severity of rotation or tilting, the inferior border of the mandible will appear distorted, and the condyles and rami will appear asymmetrical.



**Figure 4:** Incorrect lateral position tilted to the right. Note the teeth on the left are wide and poorly defined, while the teeth on the right are narrowed and severely overlapped. The inferior border of the mandible is distorted and the condyles appear asymmetrical.


To avoid imaging errors that result from incorrect lateral positioning, the midsagittal plane must be positioned perpendicular to the floor. Most panoramic x-ray machines have a head positioner and laser light beam guide, along with a mirror, to assist in determining the correct lateral head position. The pediatric patient may need additional instructions to maintain the correct position throughout the exposure.

Movement of the tube head during exposure may pique the pediatric patient's curiosity, causing the head to rotate as the eyes follow the movement of the tube head. A vertical line decal affixed to the mirror can serve as a visual aid and a focus point. An eye-catching sticker, such as those purchased from a craft store, can be adhered to the mirror in a position that aligns with the midsagittal plane. The patient can be directed to position the head so that the sticker appears at the tip of the nose and to maintain focus on this reflection throughout the exposure. Pediatric patients may find looking at themselves in the mirror entertaining and a fun way to participate in the process.<sup>9</sup>


### Superior-Inferior (Up-Down) Positioning Error

Positioning the dental arches within the superior-inferior (up-down) dimension of the focal trough can be challenging to achieve, especially with children whose smaller size reduces the distance between the shoulders and the inferior border of the chin. When the arches are positioned incorrectly in the superior-inferior direction, the image exhibits multiple distortions, including increased overlap in the premolar regions. When the arches are positioned too far up or down, the teeth will simultaneously move into a position that is too far back or too far forward, respectively, out of the focal trough.<sup>11</sup>

Positioning the arches too far superiorly produces a characteristic "frown" or flat, downward- turned appearance to the occlusal plane (**Figure 5**). The condyles flare out and off the edges of the image, and the palate appears as a widened, thick, dense radiopacity. This positioning error results in a widened appearance of the palate and obliterates the apical regions of the maxillary teeth, compromising the images of the unerupted developing dentition. As the maxillary arch tips upward, the anterior teeth tilt backward producing the same widened appearance that results from an incorrect anterior-posterior position. Positioning the arches too far inferior produces a characteristic "smile" appearance or the upward curve of the occlusal plane, with the condyles tilting inward toward the center of the image (**Figure 6**). Depending on the severity of the downward position, the vertebrae may also curve inward and appear superimposed over the condyles, and the hyoid bone may be superimposed over the mandible blocking a clear view of the erupted and unerupted mandibular teeth.



**Figure 5:** Chin positioned too far up. Note the characteristic "frown" or flat, downward-turned appearance to the occlusal plane. The widened palate obscures the view of the maxillary apices and the developing permanent dentition.



**Figure 6:** Chin positioned too far down. Note the characteristic “smile” or upward curved appearance to the occlusal plane and the hyoid bone superimposed over the mandible.

Correct positioning of the arches in the superior-inferior dimension requires that the patient stands with erect posture while tucking the chin in and down slightly, a direction that both adults and pediatric patients often find difficult to follow without specific guidance. The result is often a slumped position with the patient hunching the neck and shoulders over to place the chin on the chin rest. The vertebrae collapse causing attenuation of the x-ray beam that produces a triangular radiopacity superimposed over the mandible, and if severe, over the anterior maxillary regions as well.

Depending on the manufacturer, panoramic x-ray machines direct the operator to position the Frankfort or the canthomeatal plane parallel to the floor, or the ala-tragus line 5° down toward the floor. This is achieved by raising or lowering the chin rest so that the appropriate landmark lines up with indicators on the machine (**Figure 2**). The patient should be directed to stand in front of the panoramic x-ray machine allowing the operator to place the chin rest in a position that is slightly higher than the patient's chin. The patient is then requested to move into the overhead assembly of the machine and remain standing tall. If further adjustment is needed, it is usually to a lowered chin position. Once the patient's chin is resting on the chin rest, it is easier to move to a lower position than to a higher one. To assist with placing the chin on the chin rest while maintaining an erect posture, the pediatric patient can be directed to stand like a soldier. Most children are familiar with the straight back, chest forward tucked chin position demonstrated by military persons, and can readily mimic this stance.

#### **Further Recommendations**

Before beginning the exposure, the patient should be directed to close the lips around the bite block and to place the tongue against the palate. Leaving the lips open will create a soft tissue shadow across the teeth that that can be mistaken for caries.<sup>7</sup>

Leaving the tongue at rest during the exposure allows the radiation to easily penetrate the space of the oral cavity between the dorsal surface of the tongue and the palate,

producing a radiolucent shadow that diminishes the diagnostic quality of the radiograph (Figure 3).

"Filling in" this space with the soft tissue of the tongue can increase the quality of the image by diminishing this radiolucent shadow. When directed to place the tongue on the roof of the mouth, the pediatric patient is likely to press only the tip of the tongue against the palate. While an adult patient can usually understand what is required when directed to swallow and note the position of the tongue, a child may be directed to suck in the cheeks, which results in pushing the tongue into a position flat against the palate.<sup>7</sup>

### **Conclusion**

In addition to these guidelines for producing error-free radiographic images for the pediatric patient, panoramic machines should be evaluated periodically for accuracy. Changes may occur over time to the focal trough that interferes with the diagnostic quality of the machine.<sup>6</sup> If a decrease in image quality is noted despite following accurate patient positioning steps, the panoramic x-ray machine should be inspected and the focal trough recalibrated. The dental hygienist who is skilled in understanding general equipment operation and pediatric patient management is more likely to produce radiographic images that result in higher diagnostic yields.

### **References**

1. Public Health Service, Food and Drug Administration, American Dental Association Council on Dental Benefits Program, Council on Dental Practice, Council on Scientific Affairs. *The Selection of Patients for Dental Radiographic Examinations*. Washington, DC: United States Department of Health and Human Services; 1987 (revised 2005):1-23.
2. Pinkham PS, Casamassimo, Fields HW, McTigue DJ, Nowak AJ. *Pediatric Dentistry Infancy Through Adolescence*. 4th ed. St. Louis: Elsevier Saunders; 2005:677.
3. Thomson EM. *Panoramic radiographs and the pediatric patient*. *Dimensions of Dental Hygiene*. 2008;6(2):26-29.
4. Rushton VE, Rout J. *Panoramic Radiology*. London: Quintessence Publishing Co Ltd; 2006:17, 84-85.
5. Kullman L, Joseph B. *Quality of digital panoramic radiography in a newly established dental school*. *Swedish Dent Journal*. 2006; 30: 165-170.
6. Akarslan ZZ, Erten H, Güngör K, et al. *Common errors on panoramic radiographs taken in a dental school*. *J Contemp Dent Pract*. 2003; 4:24-34.
7. Rushton VE, Horner K, Worthington HV. *The quality of panoramic radiographs in a sample of general dental practices*. *Br Dent J*. 1999; 26: 630-633.

8. White SC, Heslop EW, Hollender LG, Mosier KM, Ruprecht A, ShROUT MK. *Parameters of radiologic care: an official report of the American Academy of Oral and Maxillofacial Radiology. Oral Surg Oral Med Oral Pathol.* 2001; 91:498-511.
9. McDonald RE, Avery DR, Dean JA. *Dentistry for the Child and Adolescent.* 8th ed. St. Louis: Elsevier Mosby; 2000:71-72.
10. Johnson ON, Thomson EM. *Essentials of Dental Radiography for Dental Assistants and Hygienists.* 8th ed. Upper Saddle River, NJ: Prentice-Hall 2007:388-397.
11. Serman N, Horrell BM, Singer, S. *High-quality panoramic radiographs. Tips and tricks. Dentistry Today.* 2003;22(1):70-73.

### 15.5.3 Setting Exposure Values to the Age Group

For more information about this topic, refer to the Appendices **15.1 Recommended X-Ray Exposure Table**.

### 15.5.4 The References Pertinent to the Potential Risks for the Pediatric Patients

#### 1) Literature

- I. ESPELID, I. MEJÀRE, K. WEERHEIJM:  
EAPD guidelines for the use of radiographs in children, P40-48. *European Journal of Pediatric Dentistry* 1/2003 Guidelines in dental radiology is designed to avoid unnecessary exposure to X-radiation and to identify individuals who may benefit from a radiographic examination. Every prescription of radiographs should be based on an evaluation of the individual patient benefit. Due to the relatively high frequency of caries among 5-year-old children, it is recommended to consider dental radiography for each child even without any visible caries or restorations. Furthermore, radiography should be considered at 8-9 years of age and then at 12-14, which is 1-2 years after the eruption of premolars and second molars. Additional bitewing controls should be based on an overall assessment of the caries activity/risk. The high-risk patient should be examined radiographically annually, while a 2-3 year interval should be considered when caries activity/risk is low. A routine survey by radiographs, except for caries, has not been shown to provide enough information to be justified considering the balance between cost (radiation and resources) and benefit.
- MICHAEL L. TAYLOR, B.SC. TOMAS KRON, PH.D., AND RICK D. FRANICH, PH.D.:

ASSESSMENT OF OUT-OF-FIELD DOSES IN RADIOTHERAPY OF BRAIN LESIONS IN CHILDREN, *Int. J. Radiation Oncology Biol. Phys.*, Vol. -, No. -, pp. 1–7, 2010 To characterize the out-of-field doses in pediatric radiotherapy and to identify simple methods by which out-of-field dose might be minimized, to reduce the risk of secondary cancers Out-of-field doses to pediatric patients can be minimized by using simple treatment

- C. THEODORAKOU, K. HORNER, K. HOWARD, A. WALKER:

Pediatric organ and effective doses in dental cone beam computed tomography Dental CBCT has been associated with higher radiation risk to the patients compared to conventional dental X-ray imaging. Several studies have investigated the radiation doses involved in dental CBCT for adults, but none has investigated pediatric doses. This study estimates the organ and effective doses to two pediatric tissue-equivalent phantoms using thermoluminescent dosimeters for three dental CBCT units and six imaging protocols. The doses to the thyroid, salivary glands and brain ranged from 0.068mSv to 1.131mSv, 0.708mSv to 2.009mSv and 0.031mSv to 1.584mSv respectively. The skin and red bone marrow have received much lower doses than the other three organs. The effective doses ranged from 0.022 mSv to 0.081 mSv. The effective doses calculated in this study were much higher than these of panoramic X-ray imaging but lower than conventional CT

- CHIYO YAMAUCHI-KAWAURA & KEISUKE FUJII & TAKAHIKO AOYAMA & SHUJI KOYAMA & MASATO YAMAUCHI:

Radiation dose evaluation in the head and neck MDCT examinations with a 6-year-old child anthropomorphic phantom, *Pediatr Radiol* (2010) 40:1206–1214 DOI 10.1007/s00247-009-1495-z

**Background:** CT examinations of the head and neck are the most commonly performed CT studies in children, raising a concern about radiation dose and their risks to children.

**Objective:** The purpose of this study was to clarify radiation dose levels for children of 6 years of age undergoing head and neck multi-detector CT (MDCT) examinations.

**Materials and methods:** Radiation doses were measured with small-sized silicon photodiode dosimeters that were implanted at various tissue and organ positions within a standard 6-year-old anthropomorphic phantom. Organ and effective

doses of brain CT were evaluated for 19 protocols in nine hospitals on various (2–320 detector rows) MDCT scanners.

**Results:** The maximum value of the mean organ dose in brain CT was 34.3 mGy for the brain. Maximum values of mean doses for the radiosensitive lens and thyroid were 32.7 mGy for a lens in brain CT and 17.2 mGy for thyroid in neck CT. The seventy-fifth percentile of effective dose distribution in brain CT was approximately the same as the diagnostic reference level (DRL) in the 2003 UK survey.

## 2) Website

For additional information on pediatric X-ray imaging, please refer to the websites below.

- <http://www.fda.gov/radiation-emittingproducts/radiationemittingproductsandprocedures/medicalimaging/ucm298899.htm>
- <http://www.imagegently.org/>

## 15.6 Abbreviations

<b>AC</b>	Alternating Current
<b>AF</b>	Auto-Focusing
<b>AMPT</b>	Adaptive layer Mode Panoramic Tomography
<b>CAN</b>	Controlled Area Network
<b>CBCT</b>	Cone-Beam Computed Tomography
<b>CMOS</b>	Complementary Metal-Oxide -Semiconductor
<b>CT</b>	Computed Tomography
<b>DAP</b>	Dose Area Product
<b>DC</b>	Direct Current
<b>DICOM</b>	Digital Imaging and Communications in Medicine
<b>EMC</b>	Electromagnetic Compatibility
<b>ENT</b>	Ear, Nose, and Throat
<b>ESD</b>	Electrostatic Discharge
<b>EUT</b>	Equipment Under Test
<b>FDD</b>	Focal spot to Detector Distance
<b>FOD</b>	Focal spot to Object Distance
<b>FOV</b>	Field of View
<b>FPD</b>	Flat Panel Detector
<b>IEC</b>	International Electro technical Commission
<b>ISO</b>	International Standards Organization
<b>LCD</b>	Liquid Crystal Display
<b>LED</b>	Light-Emitting Diode
<b>MAR</b>	Metal Artifact Reduction
<b>MPSO</b>	Multiple Portable Socket-Outlet
<b>ODD</b>	Object to Detector Distance
<b>PA</b>	Posterior / Anterior
<b>RF</b>	Radio Frequency
<b>ROI</b>	Region of Interest
<b>SID</b>	Source to Image Receptor Distance
<b>SIP</b>	Signal Input Part
<b>SOP</b>	Signal Output Part
<b>SMV</b>	Submento-Vertical
<b>STL</b>	Stereo Lithography
<b>TMJ</b>	Temporomandibular Joint
<b>UHD</b>	Ultra-High Definition





---

**Copyright by © 2018 VATECH Co., Ltd.**

---

All rights reserved.

The documentation, brand name, and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations that may be required due to technical improvement. For the most current information, contact your **VATECH** representative.

---

Tel: (+82) 1588-9510

Email: [gcs@vatech.co.kr](mailto:gcs@vatech.co.kr)

Website: [www.vatech.co.kr](http://www.vatech.co.kr)

Headquarters: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

Factory: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

---



The CE symbol grants this product compliance to the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.

---

EC Representative; Vatech Dental Manufacturing Ltd.

Suite 3, Ground Floor, Chancery House, St. Nicholas Way, Sutton, SM1 1JB  
UK

Tel: +44-0208-652-1900, +44-0208-643-7109

Fax: +44-0208-652-1909

---



# Green 16™